



QT[®] CARTS



LOGIQUIP® QUIET TREATMENT (QT) CARTS

Our exclusive line of LogiQuip **Quiet Treatment (QT) Carts** are designed for general use or specialty applications, such as Isolation, Anesthesia and Emergency, and are available in a variety of sizes. All **QT Carts** come standard with a key lock for secure storage. With soft-closing drawers and hospital grade casters, our **QT Carts** help provide a quieter environment for patient care.

QUIET FEATURES

QT® Cart Quiet Features help to provide a quieter environment for patient care with **Soft Return Drawers**, **Silent Hospital Grade Casters**, and **Quiet eLock System**.

LARGE CAPACITY DRAWERS

Large Capacity Drawers are customizable with 12 drawer front color options and three (3) drawer heights (3", 6", 9"). Every **QT Cart** comes complete with one (1) 3" and one (1) 6" divider set.

SIMPLE ACCESSORY RAIL

LogiQuip's **Simple Accessory Rail** provides the flexibility to customize your cart with several optional accessories including glove dispenser, waste basket, and sharps container set.

SECURITY OPTIONS

Several **Security Options** are available to ensure supplies are kept secure including a **Key Lock** (standard), **Breakaway Lock**, or **Quiet eLock System**.



Soft Return Drawers



Quiet eLock System



Silent Hospital Grade Casters

LIMITED 5-YEAR WARRANTY

LogiQuip products are covered by a **Limited 5-Year Warranty** from date of receipt against defect in materials or workmanship under normal use and service, excluding **Electronic Components**, **Cart Covers**, **Custom Products**, and **LogiCell Tambour Doors**. During this period, LogiQuip will, in its sole discretion, repair defects in materials or workmanship that existed when the product was received or replace the product. Visit LogiQuip.com/Warranty for more details.

SECURITY OPTIONS

SECURITY OPTIONS

To help maintain patient safety and secure supplies, **QT Carts** are available with a wide variety of lock options. From **Breakaway Locks**, which offer quick access in critical situations, to **Proximity Keyless Entry** for the most secure applications, LogiQuip has the right lock system for every situation.

AVAILABLE LOCK OPTIONS

- **Central Key Lock** - Comes standard on ALL carts.
- **Breakaway Lock** - Vertical bar for quick access during emergency situations.
- **Keyless Electronic Lock** - Touchpad for easy access.
- **HID Proximity or iClass Lock** - Touchless badge entry records access for increased security.



Central Lock
Secures All Drawers

ELOCK SPECS

- Whisper Quiet Operation
- 5,000 User Capacity
- Automatic Relock
- Key Override
- Dual Credential Capability
- Available Wi-Fi Allows IT To:
 - Remotely Add/Delete Users
 - Automatically Download Audit Trails
 - Monitor Battery Level
- Requires LockView (VLS) For Remote Connection



eLock System With Proximity Switch
(Optional)



MODEL # VLS
Lock Management Software



HID (MODEL # VPT) and
iClass Proximity Card



MODEL # VPT
HID Proximity Tag

STANDARD QT® CARTS & RAPID SHIP

STANDARD QT CARTS

Our **Standard QT Cart** comes in multiple sizes and can be equipped with a key lock, electronic lock, or proximity lock for secure storage. These general use carts are great for storing products in all patient areas and have **ABS Top Surfaces** which are light and easy to clean. All **Standard QT Carts** have a **Sliding Table Top Extension**, perfect for bedside use, with a maximum load of 17 pounds. Many **Standard QT Carts** are available immediately through our **Rapid Ship Program**.



STANDARD QT CARTS				
MODEL	SERIES	D x W x H - IN.	NUMBER OF DRAWERS	DRAWER SIZE D x W - IN.
QT35W21K3369W	35QT*	20 x 27 x 34.5	(2) 3", (1) 6", (1) 9"	15.75 x 24
QT38W24K3669W	38QT	20 X 27 X 37.5	(1) 3", (2) 6", (1) 9"	15.75 x 24
QT41W27K33669W	41QT*	20 X 27 X 40.5	(2) 3", (2) 6", (1) 9"	15.75 x 24
QT44W30K333669W	44QT*	20 X 27 X 43.5	(3) 3", (2) 6", (1) 9"	15.75 x 24
PPE QT CART				
QTP3366-WHI	QTP	21 x 20 x 34.5-46.5	(2) 3", (1) 6"	15.75 x 14

*Available through the Rapid Ship Program

PERSONALIZED EQUIPMENT CART

This **Personalized Equipment Cart** comes standard with a key lock and 4" casters. The height adjustable worktop (range of 12") is 14"D x 18"W and is the perfect size for any laptop or tablet.

RAPID SHIP PROGRAM

Our **Rapid Ship Program** allows you quick and easy access to a selection of "stocked" **QT Cart** models. These sizes are in stock and ready to ship out. All Accessories lists on pages 9-11 are also available through this program.



MODEL #
QTP3366-WHI

NARROW & DOUBLE QT® CARTS


STANDARD NARROW QT CARTS

LogiQuip **Narrow QT Carts** come standard with a key lock, 4" casters, and **ABS Top Surfaces** which are light and easy to clean. This line-up offers a **Sliding Table Top Extension**, perfect for bedside use, with a maximum load of 17 pounds. **Narrow Carts** are available in four (4) sizes with specific models available immediately through our **Rapid Ship Program**.



NARROW QT CARTS				
MODEL	SERIES	D x W x H - IN.	NUMBER OF DRAWERS	DRAWER SIZE D x W - IN.
QTN34W21K3369W	34QTN	21 x 20 x 33.5	(2) 3", (1) 6", (1) 9"	15.75 x 14
QTN37W24K3369W	37QTN	21 x 20 x 36.5	(1) 3", (2) 6", (1) 9"	15.75 x 14
QTN40W27K3369W	40QTN*	21 x 20 x 39.5	(2) 3", (2) 6", (1) 9"	15.75 x 14
QTN43W30K333669W	43QTN	21 x 20 x 42.5	(3) 3", (2) 6", (1) 9"	15.75 x 14

*Available through the Rapid Ship Program



DOUBLE QT CART				
MODEL	SERIES	D x W x H - IN.	NUMBER OF DRAWERS	DRAWER SIZE D x W - IN.
QTD34W21K3666W	34QTD	21.5 x 44 x 34	(1) 3", (3) 6"	Left: 15.75 x 24 Right 15.75 x 14

SPECIALTY QT[®] CARTS & ACCESSORY PACKAGES

EMERGENCY CARTS

Emergency QT (Crash) Carts are designed with high pressure situations in mind including the signature red color and a breakaway lock system. Using a wide assortment of accessories, these carts have the right product for any emergency situations.

Our **QT Cart Emergency Package Model QT-EMERG-PKG** includes the following accessories:

- VIV – IV Pole
- VDS – Defibrillator Shelf
- VCPR – Cardiac Board and Hooks
- V02 – Oxygen Tank Holder
- VTR – Three-Sided Rail
- AD3 – Two (2) 3" Adjustable Drawer Dividers
- AD6 – 6" Adjustable Drawer Dividers
- TOS – 15 Outlet Strip



MODEL # QT41B27K33669B
Shown with Accessories



MODEL # QT41R27B33669R
Shown with Accessories

ANESTHESIA CARTS

Anesthesia Carts store all necessary tools in a single, compact portable location with a locking system designed to keep unauthorized personnel out. The blue color makes cart identification easy.

Our **QT Cart Anesthesia Package Model QT-ANES-PKG** includes the following accessories:

- VIV – IV Pole
- P-AB – Back Bridge
- VTB-1 – Two (2) Tilt Bins, 5-Bin System
- VGD – Glove Dispenser
- VWB – Waste Basket
- AD3 – Two (2) 3" Adjustable Drawer Dividers
- AD6 – 6" Adjustable Drawer Dividers
- VTR – Three-Sided Rail
- TD-1 – Tape Dispenser

SPECIALTY QT® CARTS & ACCESSORY PACKAGES



ISOLATION CARTS

Isolation Carts are a convenient way to store PPE in an easy-access location. The bright yellow color is a constant reminder to use proper PPE and precautions for infection control when entering a patient room.

Our **QT Isolation Package Model QT-ISO-PKG** includes the following accessories:

- **VWB** – Waste Basket
- **VGD** – Glove Dispenser
- **VTR** – Three-Sided Rail
- **VUH** – Utility Hook
- **VCR** – Coat Rack

MODEL # QT44Y30K3999Y
Shown with Accessories

PEDIATRIC EMERGENCY CARTS

Pediatric Emergency Carts come with nine (9) drawers that are color code compatible with pediatric color-coding systems. A **Key Lock** and **Sliding Table Top** extension (max load of 17 pounds) are standard.

Our **QT Pediatric Model QT44W30BPEDS2** includes the following:

- Lightweight Aluminum Design
- Aluminum Rails for Simple Accessory
- Overall Dimension – 20" D x 27" W x 43.5" H
- 30" Vertical Drawer Space
- Drawer Interior Dimensions – 15.25" D x 21.75" W
- Each Drawer Has Separate Breakaway Locks
- 5" Quiet Medical Casters – One (1) Swivel, Two (2) Brake,
- One (1) Swivel/Lock







MODEL #
QT44W30BPEDS2

CUSTOM CONFIGURATIONS & COLOR OPTIONS

CUSTOM CART/CABINET CONFIGURATIONS

LogiQuip offers full customization of **QT Carts** from cabinet size, drawer size, and color.

- Pick the **Cart Body** – Standard (27" wide) or Narrow (20" wide)
- Pick the **Body Color** –  White,  Red,  Yellow,  Blue
- Pick your **Drawer Size(s)** – 3", 6", 9"
- Pick your **Drawer Color** (see **Color Options** below)
- Add your **Accessories** (see pages 9-11)



COLOR OPTIONS

LogiQuip **QT Carts** are available in twelve (12) colors. **Body Colors** in White and limited quantities of Red, Blue and Yellow are available through the **Rapid Ship Program**. All **Drawer Colors** and **Sizes** are in stock and available through **Rapid Ship**. (Actual color may vary, depending on your monitor setting or press run.)



Blue (B)



Cream (C)



Green (G)



Lavender (L)



Light Blue (LB)



Orange (O)



Pink (P)



Red (R)



Silver (S)



Taupe (T)



Yellow (Y)



White (W)

STOCK ACCESSORIES (RAPID SHIP PROGRAM)

STOCK ACCESSORIES

All Stock Accessories are available for immediate shipment through our Rapid Ship Program.



VAD-1



CPD



QT-AD3, QTN-AD3



QT-AD6, QTN-AD6



VMT-1



STOCK ACCESSORIES	
MODEL	DESCRIPTION
VAD-1	Locking Medicine Drawer
CPD	Computer Drawer
QT-AD3	3" Standard Cart – 3 short, 6 adjustable, 2-3" strips, 6 T-clips
QTN-AD3	3" Narrow Cart Set – 3 short, 6 adjustable, 2-3" strips, 6 T-clips
QT-AD6	6" Standard Cart Set – 3 short, 6 adjustable, 2-6" strips, 4 T-clips
QTN-AD6	6" Narrow Cart Set – 3 short, 6 adjustable, 2-6" strips, 4 T-clips
VMT-1	Medicine 15" x 22" x 2" with 9 compartments
VTB-1†	Large Tilt Bins (5)
VTB-2†	Medium Tilt Bins (12)

†Requires Trellis

STOCK ACCESSORIES (RAPID SHIP PROGRAM)



V02



VUH



VIV



VDS




VGB



VSC



VCPR



VBS-1†




VTS†

VAB

STOCK ACCESSORIES	
MODEL	DESCRIPTION
V02	O2 Holder and Accessory Bar
VUH	Utility Hook Rack and Accessory Bar
VIV	IV Pole
VDS	Defibrillator Shelf
VGB	Garbage Can With Lid, Basket and Accessory Bar
VWB	Waste Basket and Accessory Bar
VSC	Sharp Container Set, Basket and Accessory Bar
VAB	Accessory Bar
VCPR	Cardiac Board and Two (2) Mounting Hooks
VTS†	Top Shelf and Accessory Bar for QT Carts
VBS-1†	Small Fixed Bins and Accessory Bar Each Bin 7" x 4" x 3"

†Requires Trellis

STOCK ACCESSORIES (RAPID SHIP PROGRAM)



TD-1



VTR



VCR



P-AB



VGD-3




VFH




VGD



VHL



VSH



VMS

STOCK ACCESSORIES	
MODEL	DESCRIPTION
TD-1	Tape Dispenser
VCR	Coat Rack
P-AB	Back Bridge
VTR	Three Sided Rail
VGD-3	Triple Glove Box Holder, Accessory Bar (2)
VFH	File Holder and Accessory Bar 6" x 15" x 22"
VGD	Gloves Box Holder and Accessory Bar
VHL	Push Handle with Accessory Bar
VSH	Suction Machine and Accessory Bar
VMS	Multi-Storage Basket and Accessory Bar 7" x 13" x 5"

