



Case

CARTS



TABLE OF CONTENTS

TABLE OF CONTENTS

LOGIQUIP® CASE CARTS

FEATURES & BENEFITS 3-4

STAINLESS STEEL BOLTED CASE CARTS 5

SHELVING OPTIONS FOR BOLTED CASE CARTS 5

OVER-THE-ROAD CASE CARTS, ACCESSORIES & CASTERS 6

CASE CARTS ACCESSORIES & CASTERS 7



LOGIQUIP® CASE CARTS

Designed specifically to protect sterile supplies as they move from central processing to surgical suites, LogiQuip **Case Carts** provide a safe, durable, and efficient means of transportation. All **Stainless Steel Bolted** and **Over-the-Road Case Carts** have a full wrap-around **Bumper**, **6" Casters** for easy maneuverability, and accept height adjustable **Shelving (Shelves not included)**. LogiQuip also offers **Open Case Carts** (custom option) for a lighter option that is the ideal solution for a dedicated clean or soiled delivery system from the OR to reprocessing.

STAINLESS STEEL CONSTRUCTION

An easy-to-clean, corrosion-resistant Stainless Steel body on all **Closed** and **Over-the-Road Case Carts** provides strength and durability in the harshest environment.

6-INCH CASTERS

These durable **Castors** can withstand high temperatures and high pressure when going through an automated washer. Sealed ball bearings allow these **Castors** to quietly move through your facility.

ACCESSORIES & COMPONENTS

With a variety of bolt-on accessories and components including **Card Holders**, **Clean/Dirty Wheel Indicator**, **Handles**, and **Bumpers**, you can customize these **Case Carts** to help streamline processes.

FULLY SERVICEABLE CONSTRUCTION

Case Carts are the abused daily, which can lead to damage. With **Bolted** (vs. welded) **Construction**, LogiQuip **Case Carts** are easily serviced to reduce downtime and extend the life of your **Case Carts**.



OVER THE ROAD CASE CART



BOLTED CASE CART



OPEN CASE CART
(Available as Custom)

LIMITED 5-YEAR WARRANTY

LogiQuip products are covered by a **Limited 5-Year Warranty** from date of receipt against defect in materials or workmanship under normal use and service, excluding **Electronic Components**, **Cart Covers**, **Custom Products**, **LogiCell Tambour Doors**, **Label Holders**, and **Castors**. During this period, LogiQuip will, in its sole discretion, repair defects in materials or workmanship that existed when the product was received or replace the product. Visit [LogiQuip.com/Warranty](https://www.logiquip.com/Warranty) for more details.

LOGIQUIP® CASE CARTS



LOGIQUIP BOLTED STAINLESS STEEL CASE CARTS

BOLTED STAINLESS STEEL CASE CARTS

LogiQuip Bolted Case Carts are easily serviceable, reducing downtime and extending the life of your Case Carts. Large enough to accommodate oversized instrument sets, these Case Carts come standard with Ergonomic Handles, 270° Door Swings, 6" Casters, and Wrap Around Bumper. The Stainless Steel Cart Body is easy-to-clean, has no water retention points, and comes with a no drip "tilted" pan base. Shelving options are sold separately.

BOLTED STAINLESS STEEL CASE CARTS		
MODEL	EXTERIOR D x W x H - IN.	INTERIOR D x W x H - IN.
LQCC3630	27 x 36 x 30	25 x 30.5 x 22
LQCC3640	27 x 36 x 40	25 x 30.5 x 32
LQCC3655	27 x 36 x 55	25 x 30.5 x 47
LQCC4640	27 x 46 x 40	25 x 40.25 x 32
LQCC4655	27 x 46 x 55	25 x 40.25 x 47
LQCC5240	28 x 52 x 40	25 x 46.5 x 32
LQCC5255	28 x 52 x 55	25 x 46.5 x 47




MODEL # LQCC3640

SHELF OPTIONS FOR BOLTED CASE CARTS

There are three (3) Roll-Out Shelving options available for Bolted Stainless Steel Case Carts; Perforated Stainless Steel Shelves, Solid Stainless Steel Shelves, and Wire Case Shelves. All Shelves are sold separately and come with Adjustable Rails. Maximum weight capacity is 200 pounds for Stainless Steel Shelves and 150 pounds for the Wire Shelves.


BOLTED CASE CART SHELF OPTIONS		
MODEL	DESCRIPTION	D x W - IN.
LQCCS36PR	Roll-Out Perforated Stainless Steel Case Cart Shelf for LQCC3630, LQCC3640, and LQCC3655	23 x 26
LQCCS36SR	Roll-Out Solid Stainless Steel Case Cart Shelf for LQCC3630, LQCC3640, and LQCC3655	23 x 26
LQCCS36WR	Roll-Out Wire Case Cart Shelf for LQCC3630, LQCC3640, and LQCC3655	23 x 26
LQCCS46PR	Roll-Out Perforated Stainless Steel Case Cart Shelf for LQCC4640 and LQCC4655	23 x 36
LQCCS46SR	Roll-Out Solid Stainless Steel Case Cart Shelf for LQCC4640 and LQCC465	23 x 36
LQCCS46WR	Roll-Out Wire Case Cart Shelf for LQCC4640 and LQCC465	23 x 36
LQCCS52PR	Roll-Out Perforated Stainless Steel Case Cart Shelf for LQCC5240 and LQCC5255	23 x 42
LQCCS52SR	Roll-Out Solid Stainless Steel Case Cart Shelf for LQCC5240 and LQCC5255	23 x 42
LQCCS52WR	Roll-Out Wire Case Cart Shelf for LQCC5240 and LQCC5255	23 x 42



Roll-Out Perforated Stainless Steel Shelf



Roll-Out Solid Stainless Steel Shelf



Roll-Out Wire Shelf

OVER-THE-ROAD CASE CARTS

OVER-THE-ROAD CASE CARTS

Building on our **Standard Case Carts**, LogiQuip's **Over-The-Road Case Carts** provide a rugged, durable and safe way to transport sterile supplies when Central Processing is off-site. All welded, stainless construction, heavy-duty 6" x 2" washable casters and optional accessories including a **Tow Hitch**, help to create a clean and safe environment.

- Exterior Dimension: 30" D x 43.75" W x 57" H (8.28 sq. ft.)
- Shelf Size: 26" D x 39" W
- Capacity: Six (6) full-size containers
- Container Capacity: Available for 12, 18, and 24 containers
- Unique container retention system
- Storage Space: Top 5"H; Bottom 14.5" H
- Casters: Two (2) directional locking; Two (2) swivel with brake
- Bumper System: Full wrap around on top and bottom of cart
- Door System: Overlapping, self-sealing with 95° cart seal and double locking key
- Security Tab: Breakaway



MODEL #
OTRCC3657SS



MODEL # OTR0100-01-KT Easy slide
CLEAN/DIRTY Indicator eliminates
guesswork and errors



MODEL#
CW31-76RCSR
Rigid Caster

MODEL#
CW31-76BCSR
Swivel Caster w/Brake

OVER-THE-ROAD CASE CART CASTERS/ACCESSORIES	
MODEL	DESCRIPTION
OTR0100-01-KT	Clean/Dirty Indicator
CW31-76RCSR	Cart Caster - Rigid
CW31-76BCSR	Cart Caster - Swivel with Brake
CC-FLAG3	Flag System
OTRCC-TOW	Over-the-Road Case Cart Tow Hitch



MODEL# OTRCC-TOW
Tow Hitch System



MODEL# CC-FLAG3
Flag System

CASE CARTS ACCESSORIES & CASTERS

CLOSED CASE CART ACCESSORIES

Accessories for LogiQuip's Closed Case Carts are shown at right and below:

CLOSED CASE CART ACCESSORIES	
MODEL	DESCRIPTION
BO-0204	Pick List Card Holder - 8.5" x 11" Stainless
CC-TSE	Top Extension Drop Shelf - 16" x 24"
CC0230-01-BIO	CLEAN/DIRTY Indicator Wheel
CC-TOW	Tow Hitch System
CC-FLAG3	Flag System (shown on page 6)



MODEL# CC0230-01-BIO
CLEAN/DIRTY Indicator Wheel



MODEL# CC-TSE
Extension Drop Shelf



MODEL# BO-0204
Card Holder



CC-TOW
Tow Hitch System

CLOSED CASE CART CASTERS

Our **Standard Case Cart Casters** are 6" x 1.5" stainless steel, all swivel (one is directional locking and two are braking), and have the following:

- Sealed bearings - upper and lower
- Grease fitting - bearings and axle
- High temp wheels - poly
- Cart washable
- Optional sizes available upon request



MODEL#
CCCAST6



MODEL#
CCCAST6B



MODEL#
CCCAST6DL

CLOSED CASE CART CASTERS	
MODEL	DESCRIPTION
CCCAST6	Cart Caster - Swivel
CCCAST6B	Cart Caster - Swivel with Brake
CCCAST6DL	Cart Caster - Directional Lock



www.logiquip.com

(800) 665-3760 | customersupport@logiquip.net

8.2022