

ParStor® TOUCHLESS TRANSPORT BASKETS

TOUCHLESS WRAPPED INSTRUMENT SET TRANSPORT

Reduce the risk of **Wrapped Instrument Set** contamination with LogiQuip's **NEW ParStor Touchless Transport Basket**. This innovation allows for **Wrapped Instrument Sets** to be moved through the **Sterilization Process** without touching to help remove the risk of snagging or tearing that can lead to the contamination of **Sterilized Wrapped Instruments Sets**.

TOUCHLESS

Our **ParStor Touchless Transport Basket** allows **Wrapped Instrument Sets** to move through the **Sterilization Process** without touching the wrapped set.

ELECTROPOLISHED STAINLESS STEEL BASKETS

The **Electropolishing** process strips away embedded contaminants, burrs, micro cracks, and pits to provide a corrosion resistant **Basket** that has a smooth, bright finish.

PARSTOR & CASE CART COMPATIBLE

Our **ParStor Touchless Transport Baskets** are compatible with our **ParStor**, **ParCart**, and **Case Cart Lines**, allowing you to seamlessly move **Wrapped Instrument Sets** from **Sterile Processing** to the **OR**.

LIMITED 5-YEAR WARRANTY

LogiQuip products are covered by a **Limited 5-Year Warranty**. Visit [LogiQuip.com/Warranty](https://www.logiquip.com/Warranty) for more details.



MODEL # SPOT3N




PARSTOR® TOUCHLESS TRANSPORT BASKETS

PARSTOR TOUCHLESS TRANSPORT BASKETS AND BRACKETS

Electropolished Stainless Steel Transport Baskets and Brackets are available to retrofit an existing ParStor Cart. Up to 14 Touchless Transport Baskets can be placed in each column.

PARSTOR TRANSPORT BASKETS & BRACKETS		
MODEL	DESCRIPTION	D x W x H - IN.
SPOT3N	ParStor Touchless Transport Basket	24 x 15 x 3.5
SPOTBR	ParStor Touchless Transport Basket Brackets (Pair)	—



PARSTOR TOUCHLESS TRANSPORT BASKET KITS + ADD-ONS

LogiQuip's ParStor Movable Shelving with Touchless Transport Baskets and Brackets is available as three (3) Starter Unit Kits and three (3) Add-On Unit Kits. The height is 74.5" to the top of the Frame and 76" if you include the lip of the Top Wire Shelf. Shelving Units should be limited to rolling within the immediate area for things like cleaning. These Units are not intended to be moved between departments. All Units include four (4) Cross Bars, one (1) Wire Top Shelf, two (2) Swivel Casters, two (2) Swivel Casters with Brake and one (1) Bottom Dust Shield.

Assembly is required for all ParStor Shelving. Contact your Sales Representative to request factory or on-site assembly.

PARSTOR TRANSPORT BASKET STARTER & ADD-ON KITS		
STARTER KITS	DESCRIPTION	D x W - IN.
PS68SIMN-7	7-Tier Starter Kit, 7 Touchless Transport Baskets, 7 Sets of Transport Basket Brackets	25.5 x 21.5
PS68SIMN-8	8-Tier Starter Kit, 8 Touchless Transport Baskets, 8 Sets of Transport Basket Brackets	25.5 x 21.5
PS68SIMN-9	9-Tier Starter Kit, 9 Touchless Transport Baskets, 9 Sets of Transport Basket Brackets	25.5 x 21.5

ADD-ON KITS	DESCRIPTION	D x W - IN.
PS68SIMNA-7	Add-On Frame connects to 7-Tier Starter Kit and includes 7 Touchless Transport Baskets and 7 Sets of Transport Basket Brackets	25.5 x 19
PS68SIMNA-8	Add-On Frame connects to 8-Tier Starter Kit and includes 8 Touchless Transport Baskets and 8 Sets of Transport Basket Brackets	25.5 x 19
PS68SIMNA-9	Add-On Frame connects to 9-Tier Starter Kit and includes 9 Touchless Transport Baskets and 9 Sets of Transport Basket Brackets	25.5 x 19