



# Enclosed CARTS



# TABLE OF CONTENTS

## TABLE OF CONTENTS

|   |       |
|---|-------|
| ENCLOSED CARTS . . . . .                      | 3     |
| LOGICELL® CART FEATURES & BENEFITS . . . . .  | 4     |
| CUSTOM LOGICELL CARTS . . . . .               | 5-6   |
| DOUBLE-WIDE . . . . .                         | 5     |
| TRIPLE-WIDE . . . . .                         | 6     |
| CATHETER CARTS WITH LOGICELL . . . . .        | 7-8   |
| DOUBLE-WIDE . . . . .                         | 7     |
| TRIPLE-WIDE . . . . .                         | 8     |
| LOGICELL CARTS . . . . .                      | 9-12  |
| SINGLE-WIDE . . . . .                         | 9     |
| DOUBLE-WIDE . . . . .                         | 10-11 |
| TRIPLE-WIDE . . . . .                         | 12    |
| STAINLESS STEEL CARTS WITH LOGICELL . . . . . | 13    |
| KITS & ACCESSORIES FOR LOGICELL . . . . .     | 14-22 |
| DOUBLE-WIDE ACCESSORY KITS . . . . .          | 14-17 |
| TRIPLE-WIDE ACCESSORY KITS . . . . .          | 17-18 |
| ACCESSORIES . . . . .                         | 19-22 |
| PARSTOR® SERIES FEATURES & BENEFITS . . . . . | 23    |
| ENCLOSED PARCARTS . . . . .                   | 24    |
| QT® CARTS FEATURES & BENEFITS . . . . .       | 25    |
| STANDARD QT CARTS . . . . .                   | 26    |
| NARROW QT CARTS . . . . .                     | 27    |
| SPECIALTY & CUSTOM QT CARTS . . . . .         | 28-30 |
| QT CART ACCESSORIES . . . . .                 | 31-35 |

# LOGIQUIP® ENCLOSED CARTS

## SECURE STORAGE

LogiQuip offers several fully enclosed storage options that address safety and infection control concerns. Our LogiCell® Carts, ParCarts, and QT® Carts provide secure storage for nearly anywhere within your facility.



LogiCell Carts

ParCart

QT Carts

## LIMITED 5-YEAR WARRANTY

LogiQuip products are covered by a **Limited 5-Year Warranty** from date of receipt against defect in materials or workmanship under normal use and service, excluding **Electronic Components, Cart Covers, Custom Products, LogiCell Tambour Doors, Label Holders, and Casters**. During this period, LogiQuip will, in its sole discretion, repair defects in materials or workmanship that existed when the product was received or replace the product. Visit [LogiQuip.com/Warranty](https://LogiQuip.com/Warranty) for more details.

# LOGIQUIP® LOGICELL® CARTS

**LogiCell Carts** provide efficient enclosed storage to address safety and infection control concerns. The durable, versatile, and functional **LogiCell Carts** provide a unique **Cell Panel** along the inside walls of our **Carts** that significantly increases storage capacity. With a variety of sizes, security options, and **Accessories**, **LogiCell** is fully modular and can adapt to your needs as supplies, quantities, and processes change.

## LOGICELL PANEL

**LogiCell's** unique **Cell Panel** design allows **Tubs**, **Wire Baskets**, and **Shelves** to be placed straight in or at an angle making the **LogiCell System** incredibly efficient. In addition, **Tubs**, **Wire Baskets**, and **Shelves** can be extended 80% for maximum visibility of supplies.

## SECURITY OPTIONS

**LogiCell Carts** and **Cabinets** are available with a **Key Lock** or **eLock (Electronic Keyless Lock)** for controlled access supplies. **RFID** access is also available for hands-free access (custom orders).

## DOOR OPTIONS

**LogiCell Carts** and **Cabinets** are available with **Tambour** or **Hinged Glass Doors**. **Tambour Doors** roll into a cassette helping to keep the door clean. Our **Glass Doors** are constructed with an aluminum frame and tempered glass insert, providing visibility to the supplies inside.

## RAPID SHIP AVAILABLE

Select **LogiCell Carts** are **In-Stock NOW** and ready to ship! Contact Customer Support or your Sales Representative for more information.





# DOUBLE-WIDE CUSTOM LOGICELL® CARTS

## CUSTOM LOGICELL CARTS

**LogiCell Carts and Accessories** allow for easy customization of your **LogiCell Cart** by function, including **Cardiac**, **Laparoscopy**, **Neurology**, and more. Use a mix of **Painted Steel Shelves**, **Wire Baskets**, **Clear Tubs**, and **Catheter Slides** to configure your **LogiCell Cart** to meet your unique needs.

- Customize your **LogiCell Cart(s)** using **Shells** from pages 10-11, **Accessory Kits** from pages 14-16, and **Accessories** from pages 19-22
- The **LogiCell Carts** shown on pages 5-6 are all available with **Hinged Glass Doors**
- **Double Carts** are 28" D x 41" W x 80" H
- **Triple Carts** are 28" D x 58" W x 80" H



Cardiac Cart



Laparoscopy Cart



Neurology Cart



Ophthalmology/ENT Cart

## TRIPLE WIDE CUSTOM LOGICELL® CARTS



Arthroscopy Cart



Drape/Pack Cart



Laparoscopy Cart



Orthopedic 2 Cart

# DOUBLE-WIDE CATHETER CARTS WITH LOGICELL®

## PRE-CONFIGURED CATHETER CARTS WITH LOGICELL

- Catheter Carts come complete with **Accessories**
- **Dividers** are included with most **Tubs**, **Wire Baskets**, and **Shelves**
- Customize your **LogiCell Cart(s)** using **Shells** from pages 10-11, **Accessory Kits** from pages 14-16, and **Accessories** from pages 19-22
- **Double Carts** are 28" D x 41" W x 80" H
- **Triple Carts** are 28" D x 58" W x 80" H

| DOUBLE-WIDE CATHETER CARTS WITH LOGICELL |  |
|--|--|
| SINGLE COLUMN                            | DESCRIPTION  |
| 2651DTTGKLCAT1                           | Seven (7) 9-Hook Catheter Slides, Tambour Door                               |
| 2651DTHGKLCAT1                           | Seven (7) 9-Hook Catheter Slides, Hinged Glass Doors                         |
| 2651DTTGKLCAT2                           | Seven (7) 9-Hook Catheter Slides, Top Shelf w/5 Dividers, Tambour Door       |
| 2651DTHGKLCAT2                           | Seven (7) 9-Hook Catheter Slides, Top Shelf w/5 Dividers, Hinged Glass Doors |
| 2651DTTGKLCAT3                           | Top Shelf with 5 Dividers, Catheter Box, 22 Cells, Tambour Door              |
| 2651DTHGKLCAT3                           | Top Shelf with 5 Dividers, Catheter Box, 22 Cells, Hinged Doors              |

| DOUBLE COLUMN  | DESCRIPTION   |
|----------------|---|
| 2652DTTGKLCAT7 | Three (3) 9-hook Catheter slides, 5 Shelves w/Dividers, Tambour Door  |
| 2652DTHGKLCAT7 | Three (3) 9-hook Catheter slides, 5 Shelves w/Dividers, Hinged Glass Doors  |
| 2652DTTGKLCAT5 | Nine (9) 9-Hook Catheter Slides, Top Shelf w/2 Dividers, Tambour Door   |
| 2652DTHGKLCAT5 | Nine (9) 9-Hook Catheter Slides, Top Shelf w/2 Dividers, Hinged Glass Doors   |
| 2652DTTGKLCAT8 | 6 Shelves, one (1) 3", four (4) 5" and two (2) 8" Wire Baskets, one (1) 3", one (1) 5" and one (1) 8" Clear Plastic Tubs, Tambour Door  |
| 2652DTHGKLCAT8 | 6 Shelves, one (1) 3", four (4) 5" and two (2) 8" Wire Baskets, one (1) 3", one (1) 5" and one (1) 8" Clear Plastic Tubs, Hinged Glass Doors  |
| 2652DTTGKLCAT6 | Six (6) 9-Hook Catheter Slides, two (2) 24" Shelves - each w/2 Dividers, one (1) 3" and one (1) 5" Clear Tub w/2 Long and 2 Short Dividers, two (2) 3", two (2) 5" Wire Baskets each w/ 2 Long and 2 Short Dividers, one (1) 8" Wire Basket, Tambour Door       |
| 2652DTHGKLCAT6 | Six (6) 9-Hook Catheter Slides, two (2) 24" Shelves - each w/2 Dividers, one (1) 3" and one (1) 5" Clear Tub w/2 Long and 2 Short Dividers, two (2) 3", two (2) 5" Wire Baskets each w/ 2 Long and 2 Short Dividers, one (1) 8" Wire Basket, Hinged Glass Doors |



2651DTTGKLCAT1



2652DTHGKLCAT7



# TRIPLE-WIDE CATHETER CARTS WITH LOGICELL®

| TRIPLE-WIDE CATHETER CARTS WITH LOGICELL |  |
|--|--|
| MODEL                                    | DESCRIPTION  |
| 2653TTTGKLCAT1                           | Nine (9) 9-Hook Catheter Slides, Tambour Door  |
| 2653TTHGKLCAT1                           | Nine (9) 9-Hook Catheter Slides, Hinged Glass Doors on right side  |
| 2653TTTGKLCAT3                           | Three (3) 9-Hook Catheter Slides, 3 Metal Shelves w/2 Dividers, four (4) 3" and two (2) 5" Clear Tub w/2 Long and 2 Short Dividers, four (4) 3", two and 5" Wire Baskets each w/ 2 Long and 2 Short Dividers, one (1) 8" Wire Basket, 1 Locking Drawer, Tambour Door       |
| 2653TTHGKLCAT3                           | Three (3) 9-Hook Catheter Slides, 3 Metal Shelves w/2 Dividers, four (4) 3" and two (2) 5" Clear Tub w/2 Long and 2 Short Dividers, four (4) 3", two and 5" Wire Baskets each w/ 2 Long and 2 Short Dividers, one (1) 8" Wire Basket, 1 Locking Drawer, Hinged Glass Doors |
| 2652TTLCTGKLCAT4                         | Single Column left, Double Column right, 1 Shelf w/Dividers, two (2) 3", two (2) 5" and one (1) 8" Wire Baskets w/Dividers, one (1) 3" and one (1) 5" Clear Tub w/Dividers, Seven (7) 9-Hook Catheter Slides, Tambour Door   |
| 2652TTLCHGKLCAT4                         | Single Column left, Double Column right, 1 Shelf w/Dividers, two (2) 3", two (2) 5" and one (1) 8" Wire Baskets w/Dividers, one (1) 3" and one (1) 5" Clear Tub w/Dividers, Seven (7) 9-Hook Catheter Slides, Hinged Glass Doors on right side                             |
| 2652TTLCTGKLCAT2                         | Single Column left, Double Column right, 5 Shelves w/Dividers, two (2) 3", one (1) 5" and one (1) 8" Wire Baskets w/Dividers, two (2) 3" and one (1) 5" Clear Tub w/Dividers, Fourteen (14) 9-Hook Catheter Slides, Tambour Door   |
| 2652TTLCHGKLCAT2                         | Single Column left, Double Column right, 5 Shelves w/Dividers, two (2) 3", one (1) 5" and one (1) 8" Wire Baskets w/Dividers, two (2) 3" and one (1) 5" Clear Tub w/Dividers, Fourteen (14) 9-Hook Catheter Slides, Hinged Glass Doors on right side                       |



Model #  
2653TTTGKLCAT1

## CATHETER SLIDES

These **LogiCell Catheter Slides** are used to store catheters contained in pouches. The rack pulls out of the **Cabinet** for easy access. The front plate of each **Catheter Slide** has a 3" x 5" **Label Holder** for improved organization. **Catheter Slides** can be attached to the ceiling of a **LogiCell Cart/Cabinet**, or below a **Metal Shelf** (see page 21).



| CATHETER SLIDES |                                |
|-----------------|--------------------------------|
| MODEL           | DESCRIPTION                    |
| CAT16HD         | For 16" Shelf/Cabinet, 6 hooks |
| CAT24HD         | For 24" Shelf/Cabinet, 9 hooks |



# SINGLE-WIDE LOGICELL® CARTS

## SINGLE-WIDE, ONE COLUMN

These **Single-Wide, One Column Carts** are available in  $\frac{3}{4}$  and **Tall Heights**, are ideal for applications demanding high storage volume with limited floor space, and are sized to fit one narrow column of cell system **Accessories**. The ceiling is pre-set to hold **Catheter Slides** and is available with either **Key Lock** or **eLock**


- $\frac{3}{4}$  Height Carts 28" D x 24.5" W x 62" H
- Tall Height Carts are 28" D x 24.5" W x 80" H

| SINGLE-WIDE ONE COLUMN $\frac{3}{4}$ HEIGHT |                        |
|---|------------------------|
| MODEL                                       | DESCRIPTION            |
| 2651SQTGKL                                  | Tambour Door, Key Lock |
| 2651SQTGEL                                  | Tambour Door, eLock    |
| 2651SQHGKL                                  | Glass Door, Key Lock   |
| 2651SQHGEL                                  | Glass Door, eLock      |

| SINGLE-WIDE ONE COLUMN TALL HEIGHT |                                    |
|------------------------------------|------------------------------------|
| MODEL                              | DESCRIPTION                        |
| 2651STTGKL                         | Tambour Door, Key Lock             |
| 2651STTGEL                         | Tambour Door, eLock                |
| 2651STRGKL                         | Glass Door, Hinged Right, Key Lock |
| 2651STLGKL                         | Glass Door, Hinged Left, Key Lock  |
| 2651STRGEL                         | Glass Door, Hinged Right, eLock    |
| 2651STLGEL                         | Glass Door, Hinged Left, eLock     |



Model #  
2651SQTGKL



Model #  
2651STLGKL

## OPTIONS & ACCESSORIES

- There are 16 usable slots for **Accessories** for the  $\frac{3}{4}$  Tall Cart
- There are 22 usable slots for **Accessories** for the Tall Cart
- Carts accept **Accessories** for 24" D x 16" W columns

# DOUBLE-WIDE LOGICELL® CARTS

## DOUBLE-WIDE, ONE COLUMN

These **Double-Wide, One Column Carts** are available in  $\frac{1}{2}$ ,  $\frac{3}{4}$ , and **Tall Heights**, and are sized to fit one column of wide cell system **Accessories**. The ceiling is pre-set to hold **Catheter Slides** and is available with **Key Lock** or **eLock**.

- $\frac{1}{2}$  Height Carts are 28" D x 41" W x 42" H
- $\frac{3}{4}$  Height Carts are 28" D x 41" W x 62" H
- Tall Height Carts are 28" D x 41" W x 80" H

| DOUBLE-WIDE ONE COLUMN $\frac{1}{2}$ HEIGHT |  |
|---|--|
| MODEL                                       | DESCRIPTION                                  |
| 2651DHTGKL                                  | Tambour Door, 9 Slots/Column, Key Lock       |
| 2651DHTGEL                                  | Tambour Door, 9 Slots/Column, eLock          |
| 2651DHHGKL                                  | Hinged Glass Doors 9, Slots/Column, Key Lock |
| 2651DHHGEL                                  | Hinged Glass Doors 9, Slots/Column, eLock    |

| DOUBLE-WIDE ONE COLUMN $\frac{3}{4}$ HEIGHT |  |
|---|--|
| MODEL                                       | DESCRIPTION                            |
| 2651DQTGKL                                  | Tambour Door, 16 Slots, Key Lock       |
| 2651DQTGEL                                  | Tambour Door, 16 Slots, eLock          |
| 2651DQHGKL                                  | Hinged Glass Doors, 16 Slots, Key Lock |
| 2651DQHGEL                                  | Hinged Glass Doors, 16 Slots, eLock    |

| DOUBLE-WIDE ONE COLUMN TALL HEIGHT |  |
|------------------------------------|--|
| MODEL                              | DESCRIPTION                            |
| 2651DTTGKL                         | Tambour Door, 22 Slots, Key Lock       |
| 2651DTTGEL                         | Tambour Door, 22 Slots, eLock          |
| 2651DTHGKL                         | Hinged Glass Doors, 22 Slots, Key Lock |
| 2651DTHGEL                         | Hinged Glass Doors, 22 Slots, eLock    |



Model #  
2651DTTGEL

## OPTIONS & ACCESSORIES

- Accepts **Accessories** for 24" D x 36" W columns
- **Double Sloped Top** is available: **2650ST6D**



eLock on  
Tambour Door



eLock on  
Glass Door

# DOUBLE-WIDE LOGICELL® CARTS

## DOUBLE-WIDE, TWO COLUMNS

These **Double-Wide, Two Column Carts** are available in  $\frac{1}{2}$ ,  $\frac{3}{4}$ , and **Tall Heights**, and are sized to fit two columns of wide cell system **Accessories**. Each column is 15.75" between the cell panels. The ceiling is pre-set to hold **Catheter Slides** and is available with either **Key Lock** or **eLock**.

- **$\frac{1}{2}$  Height Carts** are 28" D x 40" W x 42" H
- **$\frac{3}{4}$  Height Carts** are 28" D x 40" W x 62" H
- **Tall Height Carts** are 28" D x 40" W x 80" H

| DOUBLE-WIDE TWO COLUMNS $\frac{1}{2}$ HEIGHT |  |
|--|--|
| MODEL  | DESCRIPTION                                  |
| 2652DHTGKL                                   | Tambour Door, 9 Slots/Column, Key Lock       |
| 2652DHTGEL                                   | Tambour Door, 9 Slots/Column, eLock          |
| 2652DHHGKL                                   | Hinged Glass Doors, 9 Slots/Column, Key Lock |
| 2652DHHGEL                                   | Hinged Glass Doors, 9 Slots/Column, eLock    |

| DOUBLE-WIDE TWO COLUMNS $\frac{3}{4}$ HEIGHT |   |
|--|---|
| MODEL  | DESCRIPTION                                   |
| 2652DQTGKL                                   | Tambour Door, 16 Slots/Column, Key Lock       |
| 2652DQTGEL                                   | Tambour Door, 16 Slots/Column, eLock          |
| 2652DQHGKL                                   | Hinged Glass Doors, 16 Slots/Column, Key Lock |
| 2652DQHGEL                                   | Hinged Glass Doors, 16 Slots/Column, eLock    |

| DOUBLE-WIDE TWO COLUMNS TALL HEIGHT |   |
|-------------------------------------|---|
| MODEL                               | DESCRIPTION                                   |
| 2652DTTGKL                          | Tambour Door, 22 Slots/Column, Key Lock       |
| 2652DTTGEL                          | Tambour Door, 22 Slots/Column, eLock          |
| 2652DTHGKL                          | Hinged Glass Doors, 22 Slots/Column, Key Lock |
| 2652DTHGEL                          | Hinged Glass Doors, 22 Slots/Column, eLock    |

## OPTIONS & ACCESSORIES

- Accepts **Accessories** for 24"D x 16" W columns
- **Double Sloped Top** is available: 2650ST6D
- **Plastic Top for Double Width  $\frac{1}{2}$  Height Cart**: 2650DHT



Model # 2652DHHGKL  
With Plastic Top Model # 2650DHT



Model # 2652DTTGEL  
Tall Height

# TRIPLE-WIDE LOGICELL® CARTS

## TRIPLE-WIDE, TWO COLUMNS

These Triple-Wide Two Column Carts are sized to fit standard cell system. Accessories are available for both Single and Double-Wide Columns. The ceiling is pre-set to hold Catheter Slides and is available with either Key Lock or eLock.

- Carts are 28" D x 58" W x 80" H

## OPTIONS & ACCESSORIES

- There are 22 usable slots for Accessories
- Accepts Accessories for 24" D x 16" W and 24" D x 36" W columns
- Triple Slope Top is available: 2650ST6T

| TRIPLE-WIDE TWO COLUMNS |   |
|-------------------------|---|
| MODEL                   | DESCRIPTION   |
| 2652TTLCTGKL            | Tambour Door, 22 Slots/Column, Small Column on Left Side, Key Lock        |
| 2652TTLCTGEL            | Tambour Door, 22 Slots/Column, Small Column on Left Side, eLock           |
| 2652TTRCTGKL            | Tambour Door, 22 Slots/Column, Small Column on Right Side, Key Lock       |
| 2652TTRCTGEL            | Tambour Door, 22 Slots/Column, Small Column on Right Side, eLock          |
| 2652TTLCHGKL            | Hinged Glass Doors, 22 Slots/Column, Small Column on Left Side, Key Lock  |
| 2652TTLCHGEL            | Hinged Glass Doors, 22 Slots/Column, Small Column on Left Side, eLock     |
| 2652TTRCHGKL            | Hinged Glass Doors, 22 Slots/Column, Small Column on Right Side, Key Lock |
| 2652TTRCHGEL            | Hinged Glass Doors, 22 Slots/Column, Small Column on Right Side, eLock    |

| TRIPLE-WIDE TWO COLUMNS SPLIT |  |
|-------------------------------|--|
| MODEL                         | DESCRIPTION  |
| 2652TTRCBSHGKL                | Hinged Glass Doors, 22 Slots/Column, Wide Column Split on Left Side, Key Lock  |
| 2652TTRCBSHGEL                | Hinged Glass Doors, 22 Slots/Column, Wide Column Split on Left, eLock          |
| 2652TTLCBSHGKL                | Hinged Glass Doors, 22 Slots/Column, Wide Column Split on Right Side, Key Lock |
| 2652TTLCBSHGEL                | Hinged Glass Doors, 22 Slots/Column, Wide Column Split on Right Side, eLock    |

| TRIPLE-WIDE THREE COLUMNS |  |
|---------------------------|--|
| MODEL                     | DESCRIPTION  |
| 2653TTTGKL                | Tambour Door, 22 Slots/Column, 3 Columns, Key Lock       |
| 2653TTTGEL                | Tambour Door, 22 Slots/Column, 3 Columns, eLock          |
| 2653TTHGKL                | Hinged Glass Doors, 22 Slots/Column, 3 Columns, Key Lock |
| 2653TTHGEL                | Hinged Glass Doors, 22 Slots/Column, 3 Columns, eLock    |



Model # 2653TTTGKL



Model # 2652TTRCTGKL  
Shown with Optional Shelves



# STAINLESS STEEL LOGICELL® CARTS

## STAINLESS STEEL SINGLE-WIDE, ONE COLUMN

These **Stainless Steel Single-Wide, One Column Carts** are sized to fit standard cell system **Accessories**, have 23 usable slots, and come standard with a **Glass Hinged Door**.

- **Carts** are 28" D x 27" W x 78" H

| SINGLE-WIDE |                                  |
|-------------|----------------------------------|
| MODEL       | DESCRIPTION                      |
| 22381SSLG   | Left Hinge Glass Door, 23 Slots  |
| 22381SSRG   | Right Hinge Glass Door, 23 Slots |



Model #  
22381SSLG

## STAINLESS STEEL DOUBLE-WIDE, TWO COLUMN

These **Stainless Steel Double-Wide, Two Column Carts** are sized to fit standard cell system **Accessories**, come standard with double **Glass Hinged Doors**, and have 23 wide slots per column.

- **Carts** are 28" D x 41" W x 78" H

| DOUBLE-WIDE |   |
|-------------|---|
| MODEL       | DESCRIPTION   |
| 24581SSHG*  | Double Glass Hinged Doors, 16 Slots/Column, 1 Column  |
| 24582SSHG   | Double Glass Hinged Doors, 16 Slots/Column, 2 Columns |

\*Can use Shelves and Catheter Slides ONLY with this Cart, not Baskets or Bins.



Model #  
24582SSHG

## STAINLESS STEEL TRIPLE-WIDE CARTS

These **Stainless Steel Triple-Wide Carts** are sized to fit standard cell system **Accessories**, have 23 usable slots per column, and come standard with triple **Glass Hinged Doors**.

- **Carts** are 28" D x 54" W x 82" H

| STAINLESS STEEL TRIPLE-WIDE |   |
|-----------------------------|---|
| MODEL                       | DESCRIPTION   |
| 24882RCSSHG                 | Triple Glass Hinged Doors, Wide Column on the left, Standard Column on the right, 23 Slots/Column |
| 24882LCSSHG                 | Triple Glass Hinged Doors, Wide Column on the right, Standard Column on the left, 23 Slots/Column |
| 24883SSHG                   | Triple Glass Hinged Doors, Three (3) Standard Columns, 23 Slots/Column (69 Slots total)           |



Model #  
24882RCSSHG

# LOGICELL® DOUBLE-WIDE CART ACCESSORY KITS

## ACCESSORY KITS FOR LOGICELL DOUBLE-WIDE, TWO COLUMN CARTS

LogiCell Accessory Kits for **Double-Wide, Two Column Carts** (21" D x 48" W x 80" H) offer multiple turnkey options thoughtfully designed by function including **Arthroscopy**, **Cardio**, **Cysto**, **Ortho**, and **Vascular**. For additional options to fit your needs, you can customize your cart with **Wire Baskets**, **Clear Tubs**, **Metal Shelves**, and **Catheter Slides** to create the ideal solution. Carts are sold separately.

## ARTHROSCOPY DOUBLE-WIDE CART ACCESSORY KIT

| ARTHROSCOPY DOUBLE-WIDE CART ACCESSORY KIT INCLUDES |   |     |                  |
|---|---|-----|------------------|
| MODEL   | DESCRIPTION                                 | QTY | D x W x H - IN.  |
| 2416SHC-2   | Steel Shelf w/Divider Slots, Frosty White   | 6   | 24 x 16 x 2.5    |
| 244DIV  | Steel Divider, Frosty White                 | 12  | 24 W x 4 H       |
| 5B  | Wire Basket, 2 Cell, White                  | 6   | 23.6 x 15.75 x 5 |
| 5BDL  | Long Wire Divider For 5B, White             | 6   | —                |
| 5BDS  | Short Wire Divider For 5B, White            | 6   | —                |
| 4TCL  | Clear Tub, 2 Cell                           | 1   | 23.6 x 15.7 x 4  |
| 4TDLCL  | Long Clear Plastic Divider For 4TCL         | 2   | —                |
| 4TDSCL  | Short Clear Plastic Divider For 4TCL        | 2   | —                |
| 8B  | Wire Basket, 3 Cell, White                  | 1   | 23.6 x 15.75 x 8 |
| 8BDL  | Long Wire Divider For 8B, White             | 1   | —                |
| 8BDS  | Short Wire Divider For B, White             | 1   | —                |
| LC-STOP   | Safety Stop for Clear Trays/Tubs (set of 2) | 1   | —                |



2652DTART  
Arthroscopy Cart

## CARDIO DOUBLE-WIDE CART ACCESSORY KIT

| CARDIO DOUBLE-WIDE CART ACCESSORY KIT INCLUDES |   |     |                  |
|--|---|-----|------------------|
| MODEL  | DESCRIPTION                                   | QTY | D x W x H - IN.  |
| 2416SHC-2                                      | Steel Shelf w/ 13 Divider Slots, Frosty White | 5   | 24 x 16 x 2.5    |
| 244DIV   | Steel Divider, Frosty White                   | 4   | 24 W x 4 H       |
| 2TCL   | Clear Tub, 1 Cell                             | 2   | 23.6 x 15.7 x 2  |
| 2TDLCL   | Long Clear Plastic Divider For 2TCL           | 4   | —                |
| 2TDSCL   | Short Clear Plastic Divider For 2TCL          | 4   | —                |
| 5B   | Wire Basket, 2 Cell, White                    | 2   | 23.6 x 15.75 x 5 |
| 5BDL   | Long Wire Divider For 5B, White               | 2   | —                |
| 5BDS   | Short Wire Divider For 5B, White              | 2   | —                |
| 4TCL   | Clear Tub, 2 Cell                             | 1   | 23.6 x 15.7 x 4  |
| 4TDLCL   | Long Clear Plastic Divider For 4TCL           | 2   | —                |
| 4TDSCL   | Short Clear Plastic Divider For 4TCL          | 2   | —                |
| 8B   | Wire Basket, 3 Cell, White                    | 2   | 23.6 x 15.75 x 8 |
| 8BDL   | Long Wire Divider For 8B, White               | 2   | —                |
| 8BDS   | Short Wire Divider For B, White               | 2   | —                |
| LC-STOP  | Safety Stop for Clear Trays/Tubs (set of 2)   | 3   | —                |



2652DCARD  
Cardio Cart

\*For Hinged Glass Doors (add HG to the end of the model number)

# LOGICELL® DOUBLE-WIDE CART ACCESSORY KITS

## CYSTO DOUBLE-WIDE CART ACCESSORY KIT

| CYSTO DOUBLE-WIDE CART ACCESSORY KIT INCLUDES |   |     |                  |  |
|---|---|-----|------------------|--|
| MODEL   | DESCRIPTION   | QTY | D x W x H - IN.  |  |
| 2416SHC-2                                     | Steel Shelf w/Divider Slots, Frosty White                     | 5   | 24 x 16 x 2.5    |  |
| 244DIV  | Steel Divider, Frosty White                                   | 12  | 24 W x 4 H       |  |
| 5B  | Wire Basket, 2 Cell, White                                    | 3   | 23.6 x 15.75 x 5 |  |
| 5BDL  | Long Wire Divider For 5B, White                               | 6   | —                |  |
| 5BDS  | Short Wire Divider For 5B, White                              | 6   | —                |  |
| 4TCL  | Clear Tub, 2 Cell   | 1   | 23.6 x 15.7 x 4  |  |
| 8B  | Wire Basket, 3 Cell, White                                    | 1   | 23.6 x 15.75 x 8 |  |
| CAT24HD                                       | Heavy-Duty, Catheter Slide Shelf, White Epoxy-Coated, 9 Hooks | 3   | 24 D             |  |
| LC-STOP                                       | Safety Stop for Clear Trays/Tubs (set of 2)                   | 1   | —                |  |



2652DCYST  
Cysto Cart

## ORTHO DOUBLE-WIDE CART ACCESSORY KIT

| ORTHO DOUBLE-WIDE CART ACCESSORY KIT INCLUDES |  |     |                     |  |
|---|--|-----|---------------------|--|
| MODEL   | DESCRIPTION                                  | QTY | D x W x H - IN.     |  |
| 2416SHC-2                                     | Steel Shelf w/13 Divider Slots, Frosty White | 4   | 24 x 16 x 2.5       |  |
| 2B  | Wire Basket, 1 Cell, White                   | 1   | 23.6 x 15.75 x 2.25 |  |
| 2BDL  | Long Wire Divider For 2B, White              | 1   | —                   |  |
| 2BDS  | Short Wire Divider For 2B, White             | 2   | —                   |  |
| 2TCL  | Clear Tub, 1 Cell                            | 1   | 23.6 x 15.7 x 2     |  |
| 2TDLCL  | Long Clear Plastic Divider For 2TCL          | 1   | —                   |  |
| 2TDSCl  | Short Clear Plastic Divider For 2TCL         | 2   | —                   |  |
| 5B  | Wire Basket, 2 Cell, White                   | 10  | 23.6 x 15.75 x 5    |  |
| 5BDL  | Long Wire Divider For 5B, White              | 10  | —                   |  |
| 5BDS  | Short Wire Divider For 5B, White             | 10  | —                   |  |
| 4TCL  | Clear Tub, 2 Cell                            | 3   | 23.6 x 15.7 x 4     |  |
| 4TDLCL  | Long Clear Plastic Divider For 4TCL          | 3   | —                   |  |
| 4TDSCl  | Short Clear Plastic Divider For 4TCL         | 3   | —                   |  |
| LC-STOP                                       | Safety Stop for Clear Trays/Tubs (set of 2)  | 4   | —                   |  |



2652DORTH  
Ortho Cart

# LOGICELL® DOUBLE-WIDE CART ACCESSORY KITS

## VASCULAR DOUBLE-WIDE CART ACCESSORY KIT

| VASCULAR DOUBLE-WIDE CART ACCESSORY KIT INCLUDES |   |     |                  |
|--|---|-----|------------------|
| MODEL  | DESCRIPTION   | QTY | D x W x H - IN.  |
| 2416SHC-2  | Steel Shelf w/Divider Slots, Frosty White                     | 3   | 24 x 16 x 2.5    |
| 244DIV   | Steel Divider, Frosty White                                   | 6   | 24 W x 4 H       |
| 5B   | Wire Basket, 2 Cell, White                                    | 2   | 23.6 x 15.75 x 5 |
| 5BDL   | Long Wire Divider For 5B, White                               | 2   | —                |
| 5BDS   | Short Wire Divider For 5B, White                              | 2   | —                |
| 4TCL   | Clear Tub, 2 Cell   | 1   | 23.6 x 15.7 x 4  |
| 4TDLCL   | Long Clear Plastic Divider For 4TCL                           | 1   | —                |
| 4TDSCL   | Short Clear Plastic Divider For 4TCL                          | 1   | —                |
| 8B   | Wire Basket, 3 Cell, White                                    | 3   | 23.6 x 15.75 x 8 |
| CAT24HD  | Heavy-Duty, Catheter Slide Shelf, White Epoxy-Coated, 9 Hooks | 3   | 24 D             |
| LC-STOP  | Safety Stop for Clear Trays/Tubs (set of 2)                   | 1   | —                |



2652DVASC  
Vascular Cart

# LOGICELL® TRIPLE-WIDE CART ACCESSORY KITS

## PRE-CONFIGURED TRIPLE-WIDE SUPPLY CARTS WITH LOGICELL

LogiCell Accessory Kits for Triple-Wide Carts (28" D x 58" W x 80" H) are available in several configurations including **Arthroscopy**, **Cardio**, **Cysto**, **Ortho**, and **Vascular**. Thoughtfully designed by function, these **Accessory Kits** make it easy to set up your **LogiCell Cart**. For additional options to fit your needs, you can customize your cart with **Wire Baskets**, **Clear Tubs**, **Metal Shelves**, and **Catheter Slides** to create the ideal solution. Carts are sold separately.

## ARTHROSCOPY TRIPLE-WIDE CART ACCESSORY KIT

| ARTHROSCOPY TRIPLE-WIDE CART ACCESSORY KIT INCLUDES |   |     |                  |
|---|---|-----|------------------|
| MODEL   | DESCRIPTION                                 | QTY | D x W x H - IN.  |
| 2416SHC-2   | Steel Shelf w/Divider Slots, Frosty White   | 8   | 24 x 16 x 2.5    |
| 244DIV  | Steel Divider, Frosty White                 | 16  | 24 W x 4 H       |
| 5B  | Wire Basket, 2 Cell, White                  | 8   | 23.6 x 15.75 x 5 |
| 4TCL  | Clear Tub, 2 Cell                           | 2   | 23.6 x 15.7 x 4  |
| 4TDLCL  | Long Clear Plastic Divider For 4TCL         | 4   | —                |
| 4TDSCL  | Short Clear Plastic Divider For 4TCL        | 4   | —                |
| 8B  | Wire Basket, 3 Cell, White                  | 3   | 23.6 x 15.75 x 8 |
| LC-STOP   | Safety Stop for Clear Trays/Tubs (set of 2) | 2   | —                |



2653TART  
Arthroscopy Cart

\*For Hinged Glass Doors (add HG to the end of the model number)



# LOGICELL® TRIPLE-WIDE CART ACCESSORY KITS

## CARDIO TRIPLE-WIDE CART ACCESSORY KIT

| CARDIO TRIPLE-WIDE CART ACCESSORY KIT INCLUDES |  |     |                  |
|--|--|-----|------------------|
| MODEL  | DESCRIPTION                                  | QTY | D x W x H - IN.  |
| 2416SHC-2                                      | Steel Shelf w/13 Divider Slots, Frosty White | 6   | 24 x 16 x 2.5    |
| 244DIV   | Steel Divider, Frosty White                  | 8   | 24 W x 4 H       |
| 2TCL   | Clear Tub, 1 Cell                            | 2   | 23.6 x 15.7 x 2  |
| 2TDLCL   | Long Clear Plastic Divider For 2TCL          | 6   | —                |
| 2TDSCL   | Short Clear Plastic Divider For 2TCL         | 6   | —                |
| 5B   | Wire Basket, 2 Cell, White                   | 5   | 23.6 x 15.75 x 5 |
| 5BDL   | Long Wire Divider For 5B, White              | 5   | —                |
| 5BDS   | Short Wire Divider For 5B, White             | 5   | —                |
| 4TCL   | Clear Tub, 2 Cell                            | 3   | 23.6 x 15.7 x 4  |
| 4TDLCL   | Long Clear Plastic Divider For 4TCL          | 3   | —                |
| 4TDSCL   | Short Clear Plastic Divider For 4TCL         | 3   | —                |
| 8B   | Wire Basket, 3 Cell, White                   | 2   | 23.6 x 15.75 x 8 |
| 8BDL   | Long Wire Divider For 8B, White              | 2   | —                |
| 8BDS   | Short Wire Divider For 8B, White             | 2   | —                |
| 8TCL   | Clear Tub, 2 Cell                            | 2   | 23.6 x 15.7 x 8  |
| 8TDLCL   | Long Wire Divider For 8TCL, White            | 4   | —                |
| 8TDSCL   | Short Wire Divider 8TCL, White               | 4   | —                |
| LC-STOP  | Safety Stop for Clear Trays/Tubs (set of 2)  | 7   | —                |



2653TCARD  
Cardio Cart

## CYSTO TRIPLE-WIDE CART ACCESSORY KIT

| CYSTO TRIPLE-WIDE CART ACCESSORY KIT INCLUDES |   |     |                     |
|---|---|-----|---------------------|
| MODEL   | DESCRIPTION   | QTY | D x W x H - IN.     |
| 2416SHC-2                                     | Steel Shelf w/Divider Slots, Frosty White                     | 7   | 24 x 16 x 2.5       |
| 244DIV  | Steel Divider, Frosty White                                   | 14  | 24 W x 4 H          |
| 2B  | Wire Basket, 1 Cell, White                                    | 6   | 23.6 x 15.75 x 2.25 |
| 2BDL  | Long Wire Divider For 2B, White                               | 6   | —                   |
| 2BDS  | Short Wire Divider For 2B, White                              | 6   | —                   |
| 2TCL  | Clear Tub, 1 Cell   | 5   | 23.6 x 15.75 x 5    |
| 2TDLCL  | Long Clear Plastic Divider For 2TCL                           | 5   | —                   |
| 2TDSCL  | Short Clear Plastic Divider For 2TCL                          | 5   | —                   |
| 5B  | Wire Basket, 2 Cell, White                                    | 3   | 23.6 x 15.75 x 5    |
| 5BDL  | Long Wire Divider For 5B, White                               | 3   | —                   |
| 5BDS  | Short Wire Divider For 5B, White                              | 3   | —                   |
| CAT24HD                                       | Heavy-Duty, Catheter Slide Shelf, White Epoxy-Coated, 9 Hooks | 3   | 24 D                |
| LC-STOP                                       | Safety Stop for Clear Trays/Tubs (set of 2)                   | 5   | —                   |



2653TCYST  
Cysto Cart

# LOGICELL® TRIPLE-WIDE CART ACCESSORY KITS

## ORTHOTRIPLE-WIDE CART ACCESSORY KIT

| ORTHOTRIPLE-WIDE CART ACCESSORY KIT INCLUDES |   |     |                     |  |
|--|---|-----|---------------------|--|
| MODEL  | DESCRIPTION                                 | QTY | D x W x H - IN.     |  |
| 2416SHC-2                                    | Steel Shelf w/Divider Slots, Frosty White   | 6   | 24 x 16 x 2.5       |  |
| 244DIV                                       | Steel Divider, Frosty White                 | 9   | 24 W x 4 H          |  |
| 2B   | Wire Basket, 1 Cell, White                  | 1   | 23.6 x 15.75 x 2.25 |  |
| 2BDL   | Long Wire Divider For 2B, White             | 1   | —                   |  |
| 2BDS   | Short Wire Divider For 2B, White            | 1   | —                   |  |
| 2TCL   | Clear Tub, 1 Cell                           | 1   | 23.6 x 15.75 x 5    |  |
| 2TDLCL                                       | Long Clear Plastic Divider For 2TCL         | 1   | —                   |  |
| 2TDSCL                                       | Short Clear Plastic Divider For 2TCL        | 1   | —                   |  |
| 5B   | Wire Basket, 2 Cell, White                  | 7   | 23.6 x 15.75 x 5    |  |
| 5BDL   | Long Wire Divider For 5B, White             | 7   | —                   |  |
| 5BDS   | Short Wire Divider For 5B, White            | 7   | —                   |  |
| 4TCL   | Clear Tub, 2 Cell                           | 4   | 23.6 x 15.7 x 4     |  |
| 4TDLCL                                       | Long Clear Plastic Divider For 4TCL         | 4   | —                   |  |
| 4TDSCL                                       | Short Clear Plastic Divider For 4TCL        | 4   | —                   |  |
| LC-STOP                                      | Safety Stop for Clear Trays/Tubs (set of 2) | 5   | —                   |  |



2653TORTH  
Ortho Cart

## VASCULAR TRIPLE-WIDE CART ACCESSORY KIT

| VASCULAR TRIPLE-WIDE CART ACCESSORY KIT INCLUDES |   |     |                     |  |
|--|---|-----|---------------------|--|
| MODEL  | DESCRIPTION                                 | QTY | D x W x H - IN.     |  |
| 2416SHC-2  | Steel Shelf w/Divider Slots, Frosty White   | 12  | 24 x 16 x 2.5       |  |
| 244DIV   | Steel Divider, Frosty White                 | 14  | 24 W x 4 H          |  |
| 2B   | Wire Basket, 1 Cell, White                  | 4   | 23.6 x 15.75 x 2.25 |  |
| 2BDL   | Long Wire Divider For 2B, White             | 8   | —                   |  |
| 2BDS   | Short Wire Divider For 2B, White            | 8   | —                   |  |
| 2TCL   | Clear Tub, 1 Cell                           | 3   | 23.6 x 15.75 x 5    |  |
| 2TDLCL   | Long Clear Plastic Divider For 2TCL         | 6   | —                   |  |
| 2TDSCL   | Short Clear Plastic Divider For 2TCL        | 6   | —                   |  |
| 5B   | Wire Basket, 2 Cell, White                  | 4   | 23.6 x 15.75 x 5    |  |
| 5BDL   | Long Wire Divider For 5B, White             | 4   | —                   |  |
| 5BDS   | Short Wire Divider For 5B, White            | 4   | —                   |  |
| 4TCL   | Clear Tub, 2 Cell                           | 3   | 23.6 x 15.7 x 4     |  |
| 4TDLCL   | Long Clear Plastic Divider For 4TCL         | 3   | —                   |  |
| 4TDSCL   | Short Clear Plastic Divider For 4TCL        | 3   | —                   |  |
| LC-STOP  | Safety Stop for Clear Trays/Tubs (set of 2) | 6   | —                   |  |



2653TVASC  
Vascular Cart

# LOGICELL® ACCESSORIES – TUBS & DIVIDERS


## CLEAR TUBS & DIVIDERS

These Polycarbonate (PC) **Clear Tubs** for **LogiCell** have an exterior dimension of 24" D x 15" W, and are used to store supplies, typically smaller or those with a concern of spills. Each **Tub** includes two (2) safety stops and can be used with standard and wide units. The **Clear Dividers** can be installed to segment the **Tub** for optimization of organization.

| TUBS*  |                                       |                 |
|--------|---------------------------------------|-----------------|
| MODEL  | DESCRIPTION                           | SLOTS USED      |
| 2TCL16 | For 16" D x 24" W Columns, With Stops | 1, 4 if slanted |
| 2TCL24 | For 24" D x 16" W Columns, With Stops | 1, 4 if slanted |
| 4TCL16 | For 24" D x 16" W Columns, With Stops | 2, 5 if slanted |
| 4TCL24 | For 16" D x 24" W Columns, With Stops | 2, 5 if slanted |
| 8TCL16 | For 24" D x 16" W Columns, With Stops | 3, 6 if slanted |
| 8TCL24 | For 16" D x 24" W Columns, With Stops | 3, 6 if slanted |

| DIVIDERS |                   |               |             |
|----------|-------------------|---------------|-------------|
| MODEL    | DESCRIPTION       | SPEC. OF TUBS | D x H - IN. |
| 2TDLCL   | For 2" Tub, Long  | Max of 8      | 22 x 2      |
| 2TDSCL   | For 2" Tub, Short | Max of 12     | 14 x 2      |
| 4TDLCL   | For 4" Tub, Long  | Max of 17     | 22 x 4      |
| 4TDSCL   | For 4" Tub, Short | Max of 27     | 14 x 4      |
| 8TDLCL   | For 8" Tub, Long  | Max of 17     | 22 x 8      |
| 8TDSCL   | For 8" Tub, Short | Max of 27     | 14 x 8      |

| TUB & DIVIDER KITS* |  |
|---------------------|--|
| MODEL               | DESCRIPTION/DIVIDERS                           |
| 2TCL16D22           | Tub for 16" D x 24" W Columns/2 Long, 2 Short  |
| 2TCL24D11           | Tub for 24" D x 16" W Columns/1 Long, 1 Short  |
| 2TCL24D12           | Tub for 24" D x 16" W Columns/1 Long, 2 Short  |
| 2TCL24D22           | Tub for 24" D x 16" W Columns/2 Long, 2 Short  |
| 2TCL24D33           | Tub for 24" D x 16" W Columns/3 Long, 3 Short  |
| 4TCL16D22           | Tub for 16" D x 24" W Columns/2 Long, 2 Short  |
| 4TCL16D33           | Tub for 16" D x 24" W Columns/3 Long, 3 Short  |
| 4TCL24D11           | Tub for 24" D x 16" W Columns/1 Long, 1 Short  |
| 4TCL24D12           | Tub for 24" D x 16" W Columns/1 Long, 2 Short  |
| 4TCL24D22           | Tub for 24" D x 16" W Columns, 2 Long, 2 Short |
| 4TCL24D33           | Tub for 24" D x 16" W Columns, 3 Long, 3 Short |
| 8TCL24D11           | Tub for 24" D x 16" W Columns/1 Long, 2 Short  |
| 8TCL24D12           | Tub for 24" D x 16" W Columns/1 Long, 2 Short  |
| 8TCL24D22           | Tub for 24" D x 16" W Columns/2 Long, 2 Short  |
| 8TCL24D33           | Tub for 24" D x 16" W Columns/3 Long, 3 Short  |



**TUB & DIVIDER KITS**  
Includes: Tub, Long Divider(s)  
and Short Divider(s)

Model # 4TCL24D11  
shown above

\*Max Weight Load – 2T Models is 20 lbs.; 4T Models is 40lbs; 8T Models 50 lbs.

# LOGICELL® ACCESSORIES – BASKETS & DIVIDERS


## WIRE BASKETS & DIVIDERS

LogiQuip's white epoxy coated steel **Wire Baskets** for **LogiCell** have an exterior dimension of 23.6" D x 15.75" W. They are used to store supplies, where the open wire structure providing air circulation is desired. Each **Wire Basket** includes two (2) built-in safety stops and can be used with standard and wide units. **Dividers** are sold individually or included in the kits to segment the **Basket** for optimization of organization.

| WIRE BASKETS* |                           |                 |
|---------------|---------------------------|-----------------|
| MODEL         | DESCRIPTION               | SLOTS USED      |
| 2B            | For 16" D x 24" W Columns | 1, 5 if slanted |
| 5B            | For 16" D x 24" W Columns | 2, 6 if slanted |
| 8B            | For 24" D x 16" W Columns | 3, 7 if slanted |

| DIVIDERS |                           |             |
|----------|---------------------------|-------------|
| MODEL    | DESCRIPTION               | D x H - IN. |
| 2BDL     | For 2" Wire Basket, Long  | 22 x 1.625  |
| 2BDS     | For 2" Wire Basket, Short | 14 x 1.625  |
| 5BDL     | For 5" Wire Basket, Long  | 22 x 4.625  |
| 5BDS     | For 5" Wire Basket, Short | 14 x 4.625  |
| 8BDL     | For 8" Wire Basket, Long  | 22 x 7.5    |
| 8BDS     | For 8" Wire Basket, Short | 14 x 7.5    |

| WIRE BASKET & DIVIDER KITS* |                 |
|-----------------------------|-----------------|
| MODEL                       | DIVIDERS        |
| 2BD11                       | 1 Long, 1 Short |
| 2BD22                       | 2 Long, 2 Short |
| 2BD33                       | 3 Long, 3 Short |
| 5BD11                       | 1 Long, 1 Short |
| 5BD22                       | 2 Long, 2 Short |
| 5BD33                       | 3 Long, 3 Short |
| 8BD11                       | 1 Long, 1 Short |
| 8BD22                       | 2 Long, 2 Short |
| 8BD33                       | 3 Long, 3 Short |



**WIRE BASKET & DIVIDER KITS**  
Includes: Basket, Long Divider(s)  
and Short Divider(s)


Model # 5BD11  
shown above

\*Max Weight Load – 2B Models is 20 lbs.; 5B Models is 40lbs; 8B Models 50 lbs.



## METAL SHELVES & DIVIDERS

These **Metal Shelves** are used to store larger or bulky items in the LogiCell® units. The shelf has a 2.5" H lip on the rear to prevent items from sliding back. The front slopes down 1" H to provide a labeling surface. The **Metal Shelf** can also be inserted in the reverse direction to retain items on the shelf in a slanted position. There are slots every 1" for attaching the **Metal Divider** to the **Shelf**, and an added channel across the rear and front that accepts **Catheter Slides**.



| METAL SHELVES WITH CHANNEL |                          |            |              |                 |                     |
|----------------------------|--------------------------|------------|--------------|-----------------|---------------------|
| MODEL                      | DESCRIPTION              | CELL SLOTS | MAX DIVIDERS | CATHETER SLIDES | D x W x H - IN.     |
| 1624SHC-2                  | For 16" D x 24" W Column | 2          | 21           | Accepts up to 3 | 16 x 24 x 2.5       |
| 2416SHC-2                  | For 24" D x 16" W Column | 2          | 13           | Accepts up to 5 | 23.5 x 15.875 x 2.5 |
| 2436SHC-2                  | For 24" D x 33" W Column | 2          | 31           | Accepts up to 7 | 23.5 x 32.75 x 2.5  |

| METAL DIVIDERS |                   |             |
|----------------|-------------------|-------------|
| MODEL          | DESCRIPTION       | W x H - IN. |
| 164DIV         | For 16" D Shelves | 14 x 4      |
| 244DIV         | For 24" D Shelves | 21.75 x 4   |

## CATHETER SLIDES

These **LogiCell Catheter Slides** are used to store catheters contained in pouches. The rack pulls out of the **Cabinet** for easy access. The front plate of each **Catheter Slide** has a 3" x 5" **Label Holder** for improved organization. **Catheter Slides** can be attached to the ceiling of a **LogiCell Cart or Cabinet**, or below a **Metal Shelf** (listed above).

| CATHETER SLIDES |                                |
|-----------------|--------------------------------|
| MODEL           | DESCRIPTION                    |
| CAT16HD         | For 16" Shelf/Cabinet, 6 hooks |
| CAT24HD         | For 24" Shelf/Cabinet, 9 hooks |




## 3-INCH LABEL HOLDERS

The LH3 3-Inch Label Holder is attached by a strap which is wrapped around a horizontal wire of the **Wire Basket** and therefore can be placed in various locations across the front of the **Wire Basket**. This 1" W x 3" H **Label Holder** is available in 8 colors (listed below) and comes in bags of 50.



Model # LH3-BLUE



(Actual color may vary, depending on your monitor setting or press run.)

| 3-INCH LABEL HOLDERS |        |
|----------------------|--------|
| MODEL                | COLOR  |
| LH3-BLUE             | Blue   |
| LH3-RED              | Red    |
| LH3-PINK             | Pink   |
| LH3-YELLOW           | Yellow |
| LH3-ORANGE           | Orange |
| LH3-GREEN            | Green  |
| LH3-PURPLE           | Purple |
| LH3                  | Opaque |

## 13.9-INCH LABEL HOLDERS

The LH13 Label Holders are 13.9" W x 2" H and are mounted by sliding onto the full width of the top edge of the **Wire Basket** or **Clear Tub**. The clear labeling area is folded over itself to create a protected space for the **Label** to slide into.



Model #  
LH132PS

| 13.9-INCH LABEL HOLDERS |       |
|-------------------------|-------|
| MODEL                   | COLOR |
| LH132PS                 | Clear |
| LH134PS                 | Clear |

## LOCKABLE DRAWERS

**Lockable Drawers** are on full-extension slides and are used to store supplies in **LogiCell Units**. Each drawer uses two (2) **Cell Slots**. A top to enclose the drawer and a surface that engages with the lock are included.

| LOCKABLE DRAWERS |                  |                     |
|------------------|------------------|---------------------|
| MODEL            | DESCRIPTION      | INT. D x W - IN.    |
| 1624LDR          | For 24" W Column | 10.25 x 18.25 x 3.8 |
| 2416LDR          | For 16" W Column | 18 x 10.4 x 3.8     |
| 2436LDR          | For 33" W Column | 18 x 27.25 x 3.8    |
| 2436SSLDR-6LD    | For 33" W Column | 21.25 x 29.5 6.5    |



Model # 2416LDR

# PARSTOR® SERIES

LogiQuip's **ParStor Series** is ideal for use in **Sterile Processing Departments (SPD)** and **Operating Rooms (OR)** to store surgical instruments, peel pouches, and sutures. Our highest-density storage system, **ParStor** is designed to be used where space is at a premium. With a variety of storage accessories that extend 100%, you will have full access to your supplies when needed. If you are looking for a mobile option, our **ParCarts** are made to move. Built on a sturdy **Dolly Platform** with **5" Casters**, large quantities of supplies can be transported over thresholds and throughout your facility with ease.

## MAXIMIZE STORAGE DENSITY

Ideal for small storage areas, our **ParStor Series** utilizes every inch of vertical storage space to maximize storage capacity. The tall, **All-Welded Frames** are outfitted with **Chrome Slides** that can be mounted **Straight** or at an **Angle**. **ParStor Frames** have **Notches** to allow **Slides** and **Accessories** to be mounted at 1" intervals.

## FLEXIBLE CONFIGURATIONS

LogiQuip's **ParStor Line** is available in multiple configurations to accommodate high-density storage in any situation. **ParStor Units** come in **Stationary** and **Movable** when you need to store supplies in place. Our **ParCarts** provide **Mobile** options when you need to move supplies from location to location. **ParCarts** come in **Open** or **Enclosed**, **Single-** and **Double-Wide**, and come in three heights.

## CUSTOMIZE WITH ACCESSORIES

**ParStor Slides** accept several **Accessories** that can be extended 100% for full access to stored supplies and adjusted without tools. **Chrome Plated Waterfall Edge Shelves** are available for instrument storage that minimizes snags and tears. In addition, our **LogiCell Wire Baskets** and **Clear Tubs** are available with and without **Dividers** to fully customize your **ParStor Unit**.

## TRANSPORT BASKET

Reduce the risk of **Wrapped Instrument Set** contamination with LogiQuip's **ParStor Transport Basket**. Move **Wrapped Instrument Sets** through the sterilization process without touching to help remove the risk of snagging or tearing that can lead to the contamination.



ParStor Mobile Cart



Open ParCart with Transport Basket



Enclosed ParCart

# ENCLOSED PARCARTS

## ENCLOSED PARCARTS


The integrated **Dolly** on **Enclosed ParCarts** creates stability, allows units to be moved between floors or departments, and is ideal for non-dedicated delivery systems. The frame is all-welded 1.5" square tubular epoxy-coated steel, enclosed with a blue epoxy-coated panel on three (3) sides and the top of the unit, and has **Vertical Slots** spaced at .4375".

**ParCarts** hold standard **ParStor 16" x 24" Baskets** and **Tubs** (see **Accessories** pages 10-11). The floor is a sturdy metal base with a **Solid Shelf** insert for infection control. Four (4) 5" quiet poly-swivel **Casters**, and **Donut Bumpers** on the corners are standard. The door(s) are 3/8" clear acrylic with (3) three hospital-grade hinges and chrome handles. **Frames** have **Vertical Slots** spaced at .4375".

All **Enclosed ParCarts** include a magnetic catch on the clear acrylic door (double models include two (2) doors). A pair of 12" **Vertical Push Handles** are included, protruding 3.125" on the left side at a height of 39-51".

**Accessories** are mounted with full-extension **Slides**. Each **Slide** is adjustable, without tools, to fit the **Slots** in the **Frame**. **Slanted Slides** can slant up to four (4) slots.

| ENCLOSED PARCARTS |                               |                  |
|-------------------|-------------------------------|------------------|
| MODEL             | DESCRIPTION                   | D x W x H - IN.  |
| 500HEN-V          | Single, 46 Slots              | 29 x 25.5 x 58.5 |
| 520HEN-V          | Double, 92 Slots              | 29 x 44 x 58.5   |
| 600HEN-V          | Single, 61 Slots              | 29 x 25.5 x 74   |
| 620HEN-V          | Double, 122 Slots             | 29 x 44 x 74     |
| 500HEN-KL         | Single, 46 Slots, Key Lock    | 29 x 25.5 x 58.5 |
| 520HEN-KL         | Double, 92, Key Lock          | 29 x 44 x 58.5   |
| 500HEN-EL         | Single, 46 Slots, eLock       | 29 x 25.5 x 58.5 |
| 520HEN-EL         | Double, 92 Slots, eLock       | 29 x 44 x 58.5   |
| 600HEN-KL         | Single, 61 Slots, Key Lock    | 29 x 25.5 x 74   |
| 620HEN-KL         | Double, 122 Slots, , Key Lock | 29 x 44 x 74     |
| 600HEN-EL         | Single, 61 Slots, eLock       | 29 x 25.5 x 74   |
| 620HEN-EL         | Double, 122 Slots, eLock      | 29 x 44 x 74     |



Model # 505EN-V

| ENCLOSED PARCART KITS |                 |  |                    |
|-----------------------|-----------------|--|--------------------|
| MODEL                 | DESCRIPTION     | INCLUDES   | D x W x H - IN.    |
| 500IV-EN-V            | IV Cart, Single | (5) Pairs of Slanted Slides with Stop, (1) 5" Basket + (2) Long Dividers, (4) 8" Baskets + (8) Long Dividers | 27.5 x 25.5 x 58.5 |
| 505EN-V               | IV Cart, Single | (6) Pairs of Slanted Slides with Stop, (1) 2" Basket, (4) 5" Baskets, (1) 8" Basket                          | 27.5 x 25.5 x 58.5 |



# LOGIQUIP® QUIET TREATMENT (QT) CARTS

Our exclusive line of LogiQuip **Quiet Treatment (QT) Carts** are designed for general use or specialty applications, such as Isolation, Anesthesia and Emergency, and are available in a variety of sizes. All **QT Carts** come standard with a **Key Lock** for secure storage. With soft-closing drawers and hospital grade casters, our **QT Carts** help provide a quieter environment for patient care.

## QUIET FEATURES

**QT® Cart Quiet Features** help to provide a quieter environment for patient care with **Soft Return Drawers**, **Silent Hospital Grade Casters**, and **Quiet eLock System**.



Soft Return Drawers

## LARGE CAPACITY DRAWERS

**Large Capacity Drawers** are customizable with 12 drawer front color options and three (3) drawer heights (3", 6", 9"). Every **QT Cart** comes complete with one (1) 3" and one (1) 6" divider set.

## SIMPLE ACCESSORY RAIL

LogiQuip's **Simple Accessory Rail** provides the flexibility to customize your cart with several optional **Accessories** including glove dispenser, wastebasket, and sharps container set.

## SECURITY OPTIONS

Several **Security Options** are available to ensure supplies are kept secure including a **Key Lock** (standard), **Proximity Lock**, **Breakaway Lock**, or **Quiet eLock System**.



Proximity Lock (left)  
Quiet eLock System (right)



# STANDARD QT® CARTS & RAPID SHIP

## STANDARD QT CARTS

Our **Standard QT Cart** comes in multiple sizes and can be equipped with a **Key Lock**, **eLock (Electronic Lock)**, or **Proximity Lock** for secure storage. These general use carts are great for storing products in all patient areas and have **ABS Top Surfaces** that are light and easy-to-clean. All **Standard QT Carts** have a **Sliding Table Top Extension**, perfect for bedside use, with a maximum load of 17 pounds. Many **Standard QT Carts** are available immediately.



| STANDARD QT CARTS |        |                     |                        |                         |
|-------------------|--------|---------------------|------------------------|-------------------------|
| MODEL             | SERIES | D x W x H - IN.     | NUMBER OF DRAWERS      | DRAWER SIZE D x W - IN. |
| QT35W21K3369W     | 35QT   | 20 x 27 x 34.5      | (2) 3", (1) 6", (1) 9" | 15.75 x 24              |
| QT38W24K3669W     | 38QT   | 20 X 27 X 37.5      | (1) 3", (2) 6", (1) 9" | 15.75 x 24              |
| QT41W27K33669W    | 41QT   | 20 X 27 X 40.5      | (2) 3", (2) 6", (1) 9" | 15.75 x 24              |
| QT44W30K333669W   | 44QT   | 20 X 27 X 43.5      | (3) 3", (2) 6", (1) 9" | 15.75 x 24              |
| PPE QT CART       |        |                     |                        |                         |
| QTP3366-WHI       | QTP    | 21 x 20 x 34.5-46.5 | (2) 3", (1) 6"         | 15.75 x 14              |

## PPE QT CART

This **PPE QT Cart** comes standard with a **Key Lock** and **4" Casters**. The **Height-Adjustable Worktop** (range of 12") is 14"D x 18"W and is the perfect size for any laptop or tablet.



MODEL # QTP3366-WHI


# NARROW & DOUBLE QT® CARTS

## NARROW QT CARTS

LogiQuip Narrow QT Carts come standard with a **Key Lock**, **4" Casters**, and **ABS Top Surfaces** which are light and easy-to-clean. This line-up offers a **Sliding Table Top Extension**, perfect for bedside use, with a maximum load of 17 pounds. **Narrow Carts** are available in four (4) sizes with specific models available immediately.



| NARROW QT CARTS  |        |                 |                        |                         |
|------------------|--------|-----------------|------------------------|-------------------------|
| MODEL            | SERIES | D x W x H - IN. | NUMBER OF DRAWERS      | DRAWER SIZE D x W - IN. |
| QTN34W21K3369W   | 34QTN  | 21 x 20 x 33.5  | (2) 3", (1) 6", (1) 9" | 15.75 x 14              |
| QTN37W24K3669W   | 37QTN  | 21 x 20 x 36.5  | (1) 3", (2) 6", (1) 9" | 15.75 x 14              |
| QTN40W27K33669W  | 40QTN  | 21 x 20 x 39.5  | (2) 3", (2) 6", (1) 9" | 15.75 x 14              |
| QTN43W30K333669W | 43QTN  | 21 x 20 x 42.5  | (3) 3", (2) 6", (1) 9" | 15.75 x 14              |



| DOUBLE QT CART |        |                 |                   |                                      |
|----------------|--------|-----------------|-------------------|--------------------------------------|
| MODEL          | SERIES | D x W x H - IN. | NUMBER OF DRAWERS | DRAWER SIZE D x W - IN.              |
| QTD34W21K3666W | 34QTD  | 21.5 x 44 x 34  | (1) 3", (3) 6"    | Left: 15.75 x 24<br>Right 15.75 x 14 |

# SPECIALTY QT® CARTS & ACCESSORY PACKAGES

## EMERGENCY CARTS

**Emergency QT (Crash) Carts** are designed with high-pressure situations in mind including the signature red color and a breakaway lock system. Using a wide assortment of **Accessories**, these carts have the right product for any emergency situations.

Our **QT Cart Emergency Package Model QT-EMERG-PKG** includes the following accessories:

- **VIV** – IV Pole
- **VDS** – Defibrillator Shelf
- **VCPR** – Cardiac Board and Hooks
- **V02** – Oxygen Tank Holder
- **VTR** – Three-Sided Rail
- **AD3** – Two (2) 3" Adjustable Drawer Dividers
- **AD6** – 6" Adjustable Drawer Dividers
- **TOS** – 15 Outlet Strip



MODEL # QT41B27K33669B  
Shown with Accessories



MODEL # QT41R27B33669R  
Shown with Accessories

## ANESTHESIA CARTS

**Anesthesia Carts** store all necessary tools in a single, compact portable location with a locking system designed to keep unauthorized personnel out. The blue color makes cart identification easy.

Our **QT Cart Anesthesia Package Model QT-ANES-PKG** includes the following accessories:

- **VIV** – IV Pole
- **P-AB** – Back Bridge
- **VTB-1** – Two (2) Tilt Bins, 5-Bin System
- **VGD** – Glove Dispenser
- **VWB** – Waste Basket
- **AD3** – Two (2) 3" Adjustable Drawer Dividers
- **AD6** – 6" Adjustable Drawer Dividers
- **VTR** – Three-Sided Rail
- **TD-1** – Tape Dispenser



# SPECIALTY QT<sup>®</sup> CARTS & ACCESSORY PACKAGES



## ISOLATION CARTS

**Isolation Carts** are a convenient way to store PPE in an easy-access location. The bright yellow color is a constant reminder to use proper PPE and precautions for infection control when entering a patient room.

Our **QT Isolation Package Model QT-ISO-PKG** includes the following accessories:

- **VWB** – Waste Basket
- **VGD** – Glove Dispenser
- **VTR** – Three-Sided Rail
- **VUH** – Utility Hook
- **VCR** – Coat Rack

MODEL # QT44Y30K3999Y  
Shown with Accessories

## PEDIATRIC EMERGENCY CARTS

**Pediatric Emergency Carts** come with nine (9) drawers that are color code compatible with pediatric color-coding systems. A **Key Lock** and **Sliding Table Top** extension (max load of 17 pounds) are standard.

Our **QT Pediatric Model QT44W30BPEDS2** includes the following:

- Lightweight Aluminum Design
- Aluminum Rails for Simple Accessory
- Overall Dimension – 20" D x 27" W x 43.5" H
- 30" Vertical Drawer Space
- Drawer Interior Dimensions – 15.25" D x 21.75" W
- Each Drawer Has Separate Breakaway Locks
- 5" Quiet Medical Casters – One (1) Swivel, Two (2) Brake,
- One (1) Swivel/Lock








MODEL #  
QT44W30BPEDS2

# QT® CART CUSTOM CONFIGURATIONS & COLOR OPTIONS

## CUSTOM CART/CABINET CONFIGURATIONS

LogiQuip offers full customization of QT Carts from **Cabinet Size**, **Drawer Size**, and **Color**.

- Pick the **Cart Body** – Standard (27" wide) or Narrow (20" wide)
- Pick the **Body Color** –  White,  Red,  Yellow,  Blue
- Pick your **Drawer Size(s)** – 3", 6", 9"
- Pick your **Drawer Color** (see **Color Options** below)
- Add your **Accessories** (see pages 9-11)



Standard QT Cart



Narrow QT Cart



Light Blue Drawer Color Option  
MODEL# QT38W24E3669LB

## COLOR OPTIONS

LogiQuip QT Carts are available in twelve (12) colors. **Body Colors** in White and limited quantities of Red, Blue and Yellow are available through the **Rapid Ship Program**. All **Drawer Colors** and **Sizes** are in stock and available. (Actual color may vary, depending on your monitor setting or press run.)



Blue (B)



Cream (C)



Green (G)



Lavender (L)



Light Blue (LB)



Orange (O)



Pink (P)



Red (R)



Silver (S)



Taupe (T)



Yellow (Y)



White (W)

# QT® CART STOCK ACCESSORIES

## STOCK ACCESSORIES

All Stock Accessories are available for immediate shipment.



VAD-1



CPD



QT-AD3, QTN-AD3



QT-AD6, QTN-AD6



| STOCK ACCESSORIES |  |
|-------------------|--|
| MODEL             | DESCRIPTION  |
| VAD-1             | Locking Medicine Drawer  |
| CPD               | Computer Drawer  |
| QT-AD3            | 3" Standard Cart – 3 short, 6 adjustable, 2-3" strips, 6 T-clips     |
| QTN-AD3           | 3" Narrow Cart Set – 3 short, 6 adjustable, 2-3" strips, 6 T-clips   |
| QT-AD6            | 6" Standard Cart Set – 3 short, 6 adjustable, 2-6" strips, 4 T-clips |
| QTN-AD6           | 6" Narrow Cart Set – 3 short, 6 adjustable, 2-6" strips, 4 T-clips   |
| VTB-1†            | Large Tilt Bins (5)  |
| VTB-2†            | Medium Tilt Bins (12)  |

†Requires Trellis

# QT® CART STOCK ACCESSORIES



VMT-1



VMT-2



VCR



P-AB



VBS-1†




VTS†

| STOCK ACCESSORIES |  |
|-------------------|--|
| MODEL             | DESCRIPTION  |
| VMT-1             | 15" x 22" x 2" with 6 Dividers + 2 Rails                 |
| VMT-2             | 15" x 22" x 4" with 6 Dividers + 2 Rails                 |
| VCR               | Coat Rack  |
| P-AB              | Back Bridge  |
| VTS†              | Top Shelf and Accessory Bar for QT Carts                 |
| VBS-1†            | Small Fixed Bins and Accessory Bar Each Bin 7" x 4" x 3" |


†Requires Trellis



## QT® CART STOCK ACCESSORIES



TD-1



VTR



VDS



VWB



VWC




VGB

| STOCK ACCESSORIES |  |
|-------------------|--|
| MODEL             | DESCRIPTION  |
| TD-1              | Tape Dispenser   |
| VTR               | Three Sided Rail   |
| VGB               | Waste Basket with Touch Top Lid, Basket, and Accessory Bar |
| VWC               | Waste Bakset with Domed Lid, Basket, and Accessory Bar     |
| VWB               | Waste Basket (No Lid) and Accessory Bar                    |
| VDS               | Defibrillator Shelf  |

## QT® CART STOCK ACCESSORIES



V02



VUH



VIV



VCPR



VAB



VSC

| STOCK ACCESSORIES |   |
|-------------------|---|
| MODEL             | DESCRIPTION                                   |
| V02               | O2 Holder and Accessory Bar                   |
| VUH               | Utility Hook Rack and Accessory Bar           |
| VIV               | IV Pole                                       |
| VAB               | Accessory Bar                                 |
| VCPR              | Cardiac Board and Two (2) Mounting Hooks      |
| VSC               | Sharp Container Set, Basket and Accessory Bar |


# QT® CART STOCK ACCESSORIES



VGD-3



VGD



VMS



VFH




VSH



VHL

| STOCK ACCESSORIES |  |
|-------------------|--|
| MODEL             | DESCRIPTION  |
| VGD-3             | Triple Glove Box Holder, Accessory Bar (2)           |
| VFH               | File Holder and Accessory Bar 6" x 15" x 22"         |
| VGD               | Gloves Box Holder and Accessory Bar                  |
| VHL               | Push Handle with Accessory Bar                       |
| VSH               | Suction Machine and Accessory Bar                    |
| VMS               | Multi-Storage Basket and Accessory Bar 7" x 13" x 5" |



[www.logiquip.com](http://www.logiquip.com)

(800) 665-3760 | [customersupport@logiquip.net](mailto:customersupport@logiquip.net)

9.2023