

SELECT QT® CARTS ARE IN-STOCK NOW!

When an urgent need arises, LogiQuip's Rapid Ship Program offers a variety of in-stock products that are ready to ship, including our QT Standard and Narrow Procedure Carts. With a variety of cart sizes and Accessories in-stock, QT Carts can be customized to meet your needs quickly. For more information contact our Customer Support Team at (800) 665-3760 or email us at customersupport@logiquip.net.

STANDARD QT® CARTS

Our Standard QT Cart can be equipped with a Key Lock, eLock (Electronic Lock), or Proximity Lock for secure storage. These general use carts are great for storing products in all patient areas and have ABS Top Surfaces that are light and easy to clean. All Standard QT Carts have a Sliding Table Top Extension, perfect for bedside use, with a maximum load of 17 pounds.

STANDARD QT CARTS				
MODEL	SERIES	D x W x H - IN.	NUMBER OF DRAWERS	DRAWER SIZE D x W - IN.
QT35W21K3369W	35QT	20 x 27 x 34.5	(2) 3", (1) 6", (1) 9"	15.75 x 24
QT41W27K33669W	41QT	20 x 27 x 40.5	(2) 3", (2) 6", (1) 9"	15.75 x 24
QT44W30K333669W	44QT	20 x 27 x 43.5	(3) 3", (2) 6", (1) 9"	15.75 x 24



35QT Series



41QT Series



44QT Series



40QTN Series

NARROW QT® CARTS

This LogiQuip Narrow QT Cart comes standard with a Key Lock, 4" Casters, a Sliding Table Top Extension (perfect for bedside use), and ABS Top Surfaces which are light and easy-to-clean. Maximum load is 17 pounds.

NARROW QT CARTS				
MODEL	SERIES	D x W x H - IN.	NUMBER OF DRAWERS	DRAWER SIZE D x W - IN.
QTN40W27K33669W	40QTN	21 x 20 x 39.5	(2) 3", (2) 6", (1) 9"	15.75 x 14

ACCESSORIES FOR QT® CARTS

All **Stock Accessories** are available for immediate shipment through our **Rapid Ship Program**. You can find a full list of more than 30 **Accessories** to help optimize your QT Cart at [LogiQuip.com](https://www.logiquip.com).

