



Warming CABINETS



WHY LOGIQUIP®?



Healthcare exclusive since 1992, LogiQuip's **Inventory Management Solutions** are designed for Healthcare Professionals looking for a better way to manage supply chain tasks. With **Expert Advice**, **Innovative Solutions**, and **360° Service** that doesn't stop after the sale, LogiQuip's sole purpose is to ensure Healthcare Professionals have the supplies they need when they need them. Partnering with you to design and implement the right solutions, we give you time back to focus on what's important, Patient Care.

EXPERT ADVICE

With an in-depth understanding of the industry and a focus on **Infection Control** and **Lean Processes**, we partner with you to ensure we design the right solutions for you.

INNOVATIVE SOLUTIONS

Designed with an eye on the horizon, LogiQuip has **Innovative Solutions** that evolve with your facility. Our flexible and efficient options help to ensure you have the supplies you need when you need them.

360° SERVICE

LogiQuip's experienced **Sales Representatives** and **Award-Winning Customer Support Team** provide unparalleled service before, during, and after the sale.



LET US DO THE HEAVY LIFTING

WITH LABOR SHORTAGES AND THE EVER-INCREASING DEMAND TO DO MORE WITH LESS, the **LogiQuip Services Team** can help to decrease that burden with turnkey solutions to help ensure you have the supplies you need when you need them.

PROJECT MANAGEMENT

From our initial meeting with you through implementation, a LogiQuip **Services Team Member** will be assigned as your **Project Manager**, working with you every step of the way. Your **Project Manager** will serve as your single point of contact during the project life cycle. Together we will develop and implement project milestones until the work is complete. Contact your **Sales Representative** or call **(800) 665-3760** for more information.

ON-SITE ANALYSIS

Working with your LogiQuip **Representative**, an on-site analysis will be done to ensure that the **Project Manager** has a complete understanding of your current environment. During our on-site analysis we will review your current state, identify gaps, and align on the goals of the project.

PLANNING

Using a mix of innovative products, Lean principles, tried and true methods, and out of the box thinking, LogiQuip will provide a holistic plan with turnkey solutions to minimize downtime while decreasing your storage footprint and increasing PAR levels.

INSTALLATION & CONVERSION

LogiQuip is dedicated to providing on-time delivery, efficient installation, and supply conversion to minimize down time and time on site to get you up and running as quickly as possible.

ONGOING SUPPORT

With the project complete, we want to help you sustain the changes that have been made. With **Annual Inventory Reviews**, **Preventative Maintenance Plans**, and **Field Service** for in warranty and out of warranty repairs, we have you covered to keep things running smoothly.

LOGIQUIP® WARMING CABINETS

Increase patient satisfaction, provide better patient outcomes, and even reduce the amount of time a patient spends in the hospital with LogiQuip's **Fluid** and **Blanket Warmers**. With a variety of temperature-controlled **Warming Cabinets** to choose from, you can be sure that your **Blankets** and **IV** and **Irrigation Fluids** are stored at the optimal temperature at all times.

DATA LOGGING

Each **Fluid** and **Warming Cabinet** has an integrated **Data Logger** that anatomically record the interior temperature of the cabinet every 30 minutes. This data is stored for up to two (2) years and can be easily downloaded by simply inserting a USB drive. No additional software is needed.

TEMPERATURE CONTROL

Ideal for warming **IV Fluids**, **Solutions**, and **Blankets**, these fully insulated **Cabinets** hold their temperature to ± 3 degrees. With easy-to-use digital **Temperature Controls** that can be locked and large displays that show the actual temperature and set points at all times, you can be sure that the supplies you need are at the optimal temperature.

STANDARD FEATURES

Stainless Steel Construction, **Adjustable Feet**, **Audible** and **Visual Over-Temp Alarm**, and **Door Locks**, come standard with every **Warmer**.

OPTIONAL FEATURES

With an impressive list of optional features including **Mobile Bases** and **Stands**, **Seismic Brackets**, **Stainless Steel Roll-Out Baskets**, and **Pass Through Chambers** and **Doors**. We have the everything you need to create the right solution.

LIMITED 5-YEAR WARRANTY

LogiQuip products are covered by a **Limited 5-Year Warranty** from date of receipt against defect in materials or workmanship under normal use and service, excluding **Electronic Components**, **Cart Covers**, **Custom Products**, **LogiCell Tambour Doors**, **Label Holders**, and **Casters**. During this period, LogiQuip will, in its sole discretion, repair defects in materials or workmanship that existed when the product was received or replace the product. Visit LogiQuip.com/Warranty for more details.

DATE / TIME	TEMPERATURE	°C/°F
3-10-2019 / 23:12	106	F
3-10-2019 / 23:43	106	F
3-11-2019 / 0:12	106	F
3-10-2019 / 0:43	106	F
3-10-2019 / 1:12	106	F

Sample Report for WB2630TS Single, S/N:
062112-D01, 031115



MODEL #
W2024MDG



LOGIQUIP® WARMING CABINETS

BLANKET & FLUID WARMING CABINETS

LogiQuip's Blanket and Fluid Warmers help to provide a better experience for patients when faced with a hospital stay. Our **Data Logging Blanket and Fluid Warming Cabinets** are our most user-friendly temperature recording warmer with a simple plug and play data download that requires no additional software.

STANDARD FEATURES

- 1 Integrated Data Logger with Removable USB
- 2 USB Port
- 3 Fully Insulated Stainless Steel Cabinet
- 4 Digital Temperature Display(s)
- 5 Temperature Control(s)
- 6 Individual Chamber Power Control(s)
- 7 Keyed Temperature Lock Out
- 8 Adjustable Feet
- 9 4-Inch Base
 - Audible and Visual Over-Temp Alarm



OPTIONAL FEATURES

- A Full View Glass Door
- B Left Hand Door Hinge
- C Stainless Steel Roll Out Basket(s) with Cool Touch Handles
 - Mobile Base or Stand
 - Seismic Brackets
 - 220/240V, 50/60Hz
 - Pass Through Chambers and Doors



MODEL #
W2630TDG-LH



Scan QR Code to see
product details.

WARMING CABINETS STANDARD MODELS


WARMING CABINET MODELS

WARMING CABINET BASE MODELS					
STAINLESS STEEL DOOR	GLASS DOOR	D x W x H - IN.	CHAMBERS / SPLIT (SEE CABINET SPLITS BELOW)	# OF SHELVES	BASE
W1518XSS	W1518XSSG	17.5 x 18 x 22	Short, Single Chamber	0	—
W2024XSS	W2024XSSG	20.5 x 24 x 24.5	Short, Single Chamber	0	—
W2630XSS	W2630XSSG	26.5 x 30 x 24.5	Short, Single Chamber	1	—
W2630SS	W2630SSG	26.5 x 30 x 36	Short, Single Chamber	2	2"
W2024MS	W2024MSG	20.5 x 24 x 64.75	Medium Height, Single Chamber	2	4"
W2030MS	W2030MSG	20.5 x 30 x 64.75	Medium Height, Single Chamber	2	4"
WB2630TS	WB2630TSG	20.5 x 30 x 74.5	Tall, Single Chamber	2	4"
W2024MD	W2024MDG	20.5 x 24 x 64.75	Medium, Dual Chamber	3	4"
W2630MD	W2630MDG	20.5 x 30 x 64.75	Medium, Dual Chamber	3	4"
W2624TD	W2624TDG	26.5 x 24 x 74.5	Tall, Dual Chamber / $\frac{1}{3}$ $\frac{2}{3}$	3	4"
W2630TD	W2630TDG	26.5 x 30 x 74.5	Tall, Dual Chamber / $\frac{1}{3}$ $\frac{2}{3}$	3	4"
W2630MSTD	W2630MSTDG	26.5 x 30 x 74.5	Middle Split, Tall, Dual Chamber / $\frac{1}{2}$ $\frac{1}{2}$	4	4"
W2630TT	W2630TTG	26.5 x 30 x 78.5	Tall, Triple Chamber / $\frac{1}{3}$ $\frac{1}{3}$ $\frac{1}{3}$	1	4"


NOTE: Due to the vast number of sizes we have available, not all our models are shown. Please let us know your specific dimensions to get the model that will best suit your needs. If we don't currently have the size you need, we can custom build one for you!

SINGLE, DUAL, TRIPLE CHAMBERS


LogiQuip Warming Cabinets come with multiple cabinet configurations as shown below.



Tall, Dual Chamber



Middle Split, Tall,
Dual Chamber



Tall, Triple Chamber

OPTIONAL ADD-ONS & ACCESSORIES

WARMING CABINET AD-ONS

Warming Cabinets have a number of **Add-Ons** such as **Left Hinges**, **2 or 6-inch Bases**, **Sloped Tops** and more, that allow for customization. Simply take the base model number of the **Warming Cabinet** you want from page 6 and add the **Add-On** designation to the back of the model number. See examples below. For more information contact your **Sales Representative**.


WARMING CABINET ADD-ONS				
DOOR TYPE	DOOR HINGE	HEADER CONTROL	BASE STYLE	ADDITIONAL OPTIONS
Choose from Stainless Steel or Glass. (See chart on page 6)	Standard: Right Hinge Optional: LH = Left Hinge	Standard: Digital [†] Optional: TS = Touch Screen [†] (Single and Double only)	Standard: 4" Base [†] Optional: 2B = 2" Base [†] 6B = 6" Base [†] NB = No Base [†] M = Mobile Base MS = Mobile Stand ^{††}	ST = Slope Top [†] LDL = Lower Door Lock [†] 220 = 220 Volt Power Supply

[†] Not Available for W1518XSS or W1518XSSG

^{††} Not Available for W2024MS, W2024MSG, WB2630TS, WB2630TSG and all Dual and Triple Chambers

EXAMPLE MODELS WITH ADD-ONS	
MODEL	DESCRIPTION
W1518XSS-LH-MS	Short, Single Chamber with Stainless Steel Door; Left Hinge and Mobile Stand
W2630XSSG-TS-2B	Short, Single Chamber with Glass Door, Touch Screen and 2" Base
W2624TDG-ST-LDL	Tall, Dual Chamber / $\frac{1}{3}$ $\frac{2}{3}$ with Glass Door, Slope Top and Lower Door Lock

Mobile Base (M)



MODEL # MB2424

Mobile Stand (MS)




MODEL # MS2430

Basket Accessory



Shelf Accessory





www.logiquip.com

(800) 665-3760 | customerxxperience@logiquip.net

1.2024