

# Dri-Stor MAINTENANCE PLANS<sup>††</sup>

LogiQuip is committed to customer satisfaction throughout the life of your products. With this in mind, we are pleased to offer Maintenance Plans to keep your new LogiQuip Dri-Stor Drying Cabinet running like new. These hassle-free plans take the maintenance of your Dri-Stor Drying Cabinets off of your to-do list. Contact your LogiQuip Sales Representative for a quote.

## STANDARD

### WARRANTY

LogiQuip Dri-Stor Cabinets are covered by a **5-Year Limited Warranty** from date of receipt against defect in materials or workmanship under normal use and service, excluding electronic components. During this period, LogiQuip will, in its sole discretion, repair defects in materials or workmanship that existed when the product was received or replace the product. If the product is replaced, the new product will be warranted for the remaining term of the original product warranty period.

**1-Year Warranty** on electronic components and sub-assemblies.

## 2-YEAR

### MAINTENANCE + STANDARD WARRANTY

#### MODEL # LQDS2MP 5-YEAR LIMITED WARRANTY + 24 MONTHS OF MAINTENANCE

- Maintenance Visits Every 6 Months
- Replacement of Consumable Parts (40000CPKIT-1 or 40000CPKIT-2) Every 6 Months
- Replacement of Scope Cabinet HEPA Filter (40000MFF) Every 2 Years
- Inspect Tubing
- Check Operation of Fans and Pumps
- Check Doors and Locks
- Check Levelness of Cabinets

## 4-YEAR

### MAINTENANCE + STANDARD WARRANTY

#### MODEL # LQDS4MP 5-YEAR LIMITED WARRANTY + 48 MONTHS OF MAINTENANCE

- Maintenance Visits Every 6 Months
- Replacement of Consumable Parts (40000CPKIT-1 or 40000CPKIT-2) Every 6 Months
- Replacement of Scope Cabinet HEPA Filter (40000MFF) Every 2 Years
- Inspect Tubing
- Check Operation of Fans and Pumps
- Check Doors and Locks
- Check Levelness of Cabinets

## MAINTENANCE SCHEDULE

TYPE OF SERVICE	MODEL NUMBER	6 MONTHS	12 MONTHS
Cabinet Pre-Filter	40000EF	Replace	Replace
Channel Purge HEPA Filter Cannister (if equipped)	40000CPHF	Replace	Replace
HEPA Pre- and Post-Filters (including O-Rings)	40000PPFK	Replace	Replace
Scope Cabinet HEPA Filter <sup>*</sup>	40000MFF	Replace Every 2 years	
Tubing <sup>†</sup>	Model # Varies	Inspect	Inspect

<sup>\*</sup>Can be ordered separately

<sup>†</sup>Tubing replaced is an additional cost

<sup>††</sup>Does not apply to KP Kaiser Permanente



MODEL #  
40000MFF

# Dri-Stor MAINTENANCE KITS

Keep your Dri-Stor Drying Cabinet running like new with our Filter Maintenance Kits. Choose from three (3) Kits (below) for a one-year supply of Filters and O-Rings which need to be replaced every six (6) months.

## Kit 1

### 1-YEAR SUPPLY

Replace Filters and O-Rings  
Every 6 Months

**FOR CHANNEL PURGE  
CABINETS WITH  
10 SCOPES/PROBES OR LESS**

#### Model 40000CPKIT-1

#### INCLUDES:

**40000CPHF** – 2 Channel Purge  
Hepa Filter Canisters

**40000HEPAKIT** – 2 HEPA Cabinet  
External Filters, 2 HEPA Pre-Filters,  
2 HEPA Post-Filters, 2 O-Rings



## Kit 2

### 1-YEAR SUPPLY

Replace Filters and O-Rings  
Every 6 Months

**FOR CHANNEL PURGE  
CABINETS WITH  
11-20 SCOPES**

#### Model 40000CPKIT-2

#### INCLUDES:

**40000CPHF** – 4 Channel Purge  
Hepa Filter Canisters

**40000HEPAKIT** – 2 HEPA Cabinet  
External Filters, 2 HEPA Pre-Filters,  
2 HEPA Post-Filters, 2 O-Rings



## HEPA Kit

### 1-YEAR SUPPLY

Replace Filters and O-Rings  
Every 6 Months

**FOR HEPA FILTER MODELS  
WITHOUT CHANNEL PURGE**

#### Model 40000HEPAKIT

#### INCLUDES:

2 HEPA Cabinet External Filters

2 HEPA Pre-Filters

2 HEPA Post-Filters

2 O-Rings

