



# QT<sup>®</sup> CARTS



# WHY LOGIQUIP®?



Healthcare exclusive since 1992, LogiQuip's **Inventory Management Solutions** are designed for Healthcare Professionals looking for a better way to manage supply chain tasks. With **Expert Advice**, **Innovative Solutions**, and **360° Service** that doesn't stop after the sale, LogiQuip's sole purpose is to ensure Healthcare Professionals have the supplies they need when they need them. Partnering with you to design and implement the right solutions, we give you time back to focus on what's important, Patient Care.

## **EXPERT ADVICE**

With an in-depth understanding of the industry and a focus on **Infection Control** and **Lean Processes**, we partner with you to ensure we design the right solutions for you.

## **INNOVATIVE SOLUTIONS**

Designed with an eye on the horizon, LogiQuip has **Innovative Solutions** that evolve with your facility. Our flexible and efficient options help to ensure you have the supplies you need when you need them.

## **360° SERVICE**

LogiQuip's experienced **Sales Representatives** and **Award-Winning Customer Support Team** provide unparalleled service before, during, and after the sale.



WITH LABOR SHORTAGES AND THE EVER-INCREASING DEMAND TO DO MORE WITH LESS, the **LogiQuip Services Team** can help to decrease that burden with turnkey solutions to help ensure you have the supplies you need when you need them.

## PROJECT MANAGEMENT

From our initial meeting with you through implementation, a LogiQuip **Services Team Member** will be assigned as your **Project Manager**, working with you every step of the way. Your **Project Manager** will serve as your single point of contact during the project life cycle. Together we will develop and implement project milestones until the work is complete. Contact your **Sales Representative** or call **(800) 665-3760** for more information.

## ON-SITE ANALYSIS

Working with your LogiQuip **Representative**, an on-site analysis will be done to ensure that the **Project Manager** has a complete understanding of your current environment. During our on-site analysis we will review your current state, identify gaps, and align on the goals of the project.

## PLANNING

Using a mix of innovative products, Lean principles, tried and true methods, and out of the box thinking, LogiQuip will provide a holistic plan with turnkey solutions to minimize downtime while decreasing your storage footprint and increasing PAR levels.

## INSTALLATION & CONVERSION

LogiQuip is dedicated to providing on-time delivery, efficient installation, and supply conversion to minimize down time and time on site to get you up and running as quickly as possible.

## ONGOING SUPPORT

With the project complete, we want to help you sustain the changes that have been made. With **Annual Inventory Reviews**, **Preventative Maintenance Plans**, and **Field Service** for in warranty and out of warranty repairs, we have you covered to keep things running smoothly.

# TABLE OF CONTENTS

## TABLE OF CONTENTS

WHY LOGIQUIP®? . . . . .	2
LOGIQUIP SERVICES TEAM . . . . .	3
LOGIQUIP QT® (QUIET TREATMENT) CARTS	
FEATURES & BENEFITS . . . . .	5
STANDARD QT CARTS. . . . .	6
NARROW QT CARTS . . . . .	7
SPECIALTY QT CARTS & ACCESSORY PACKAGES . . . . .	8-9
SECURITY OPTIONS . . . . .	10
CUSTOM CONFIGURATIONS & COLOR OPTIONS . . . . .	11
STOCK ACCESSORIES (RAPID SHIP PROGRAM) . . . . .	12-15





# LOGIQUIP® QUIET TREATMENT (QT) CARTS

Our exclusive line of LogiQuip **Quiet Treatment (QT) Carts** are designed for general use or specialty applications, such as Isolation, Anesthesia and Emergency, and are available in a variety of sizes. All **QT Carts** come standard with a **Key Lock** for secure storage. With soft-closing drawers and hospital grade casters, our **QT Carts** help provide a quieter environment for patient care.

## QUIET FEATURES

**QT® Cart Quiet Features** help to provide a quieter environment for patient care with **Soft Return Drawers**, **Silent Hospital Grade Casters**, and **Quiet eLock System**.

## LARGE CAPACITY DRAWERS

**Large Capacity Drawers** are customizable with 12 drawer front color options and three (3) drawer heights (3", 6", 9"). Every **QT Cart** comes complete with one (1) 3" and one (1) 6" divider set.

## SIMPLE ACCESSORY RAIL

LogiQuip's **Simple Accessory Rail** provides the flexibility to customize your cart with several optional **Accessories** including glove dispenser, wastebasket, and sharps container set.

## SECURITY OPTIONS

Several **Security Options** are available to ensure supplies are kept secure including a **Key Lock** (standard), **Breakaway Lock**, or **Quiet eLock System**.



## LIMITED 5-YEAR WARRANTY

LogiQuip products are covered by a **Limited 5-Year Warranty** from date of receipt against defect in materials or workmanship under normal use and service, excluding **Electronic Components**, **Cart Covers**, **Custom Products**, **LogiCell Tambour Doors**, **Label Holders**, and **Casters**. During this period, LogiQuip will, in its sole discretion, repair defects in materials or workmanship that existed when the product was received or replace the product. Visit [LogiQuip.com/Warranty](https://www.logiquip.com/Warranty) for more details.

# STANDARD QT® CARTS & RAPID SHIP

## STANDARD QT CARTS

Our Standard QT Cart comes in multiple sizes and can be equipped with a **Key Lock**, **eLock (Electronic Lock)**, or **Proximity Lock** for secure storage. These general use carts are great for storing products in all patient areas and have **ABS Top Surfaces** that are light and easy-to-clean. All Standard QT Carts have a **Sliding Table Top Extension**, perfect for bedside use, with a maximum load of 17 pounds. Many Standard QT Carts are available immediately through our **Rapid Ship Program**.



STANDARD QT CARTS				
MODEL	SERIES	D x W x H - IN.	NUMBER OF DRAWERS	DRAWER SIZE D x W - IN.
QT35W21K3369W	35QT*	20 x 27 x 34.5	(2) 3", (1) 6", (1) 9"	15.75 x 24
QT38W24K3669W	38QT	20 X 27 X 37.5	(1) 3", (2) 6", (1) 9"	15.75 x 24
QT41W27K33669W	41QT*	20 X 27 X 40.5	(2) 3", (2) 6", (1) 9"	15.75 x 24
QT44W30K333669W	44QT*	20 X 27 X 43.5	(3) 3", (2) 6", (1) 9"	15.75 x 24
PPE QT CART				
QTP3366-WHI	QTP	21 x 20 x 34.5-46.5	(2) 3", (1) 6"	15.75 x 14

\*Available through the Rapid Ship Program

## PPE QT CART

This PPE QT Cart comes standard with a **Key Lock** and **4" Casters**. The **Height-Adjustable Worktop** (range of 12") is 14"D x 18"W and is the perfect size for any laptop or tablet.

## RAPID SHIP PROGRAM

Our **Rapid Ship Program** allows you quick and easy access to a selection of "stocked" QT Cart models. Select sizes are in stock and ready to ship out. All **Accessories** listed on pages 13-15 are also available through this program.



MODEL # QTP3366-WHI

# NARROW & DOUBLE QT® CARTS


## NARROW QT CARTS

LogiQuip Narrow QT Carts come standard with a **Key Lock**, **4" Casters**, and **ABS Top Surfaces** which are light and easy-to-clean. This line-up offers a **Sliding Table Top Extension**, perfect for bedside use, with a maximum load of 17 pounds. **Narrow Carts** are available in four (4) sizes with specific models available immediately through our **Rapid Ship Program**.



NARROW QT CARTS				
MODEL	SERIES	D x W x H - IN.	NUMBER OF DRAWERS	DRAWER SIZE D x W - IN.
QTN34W21K3369W	34QTN	21 x 20 x 33.5	(2) 3", (1) 6", (1) 9"	15.75 x 14
QTN37W24K33669W	37QTN	21 x 20 x 36.5	(1) 3", (2) 6", (1) 9"	15.75 x 14
QTN40W27K33669W	40QTN*	21 x 20 x 39.5	(2) 3", (2) 6", (1) 9"	15.75 x 14
QTN43W30K333669W	43QTN	21 x 20 x 42.5	(3) 3", (2) 6", (1) 9"	15.75 x 14

\*Available through the Rapid Ship Program



DOUBLE QT CART				
MODEL	SERIES	D x W x H - IN.	NUMBER OF DRAWERS	DRAWER SIZE D x W - IN.
QTD34W21K3666W	34QTD	21.5 x 44 x 34	8 Soft Close Drawers* Left: (Wide) one (1) 3" and three (3) 6" Right: (Narrow) one (1) 3" and three (3) 6"	Left (Wide): 15.75 x 24 Right (Narrow): 15.75 x 14

\*Includes two (2) 3" and two (2) 6" Drawer Divider Sets.

# SPECIALTY QT® CARTS & ACCESSORY PACKAGES

## EMERGENCY CARTS

Emergency QT (Crash) Carts are designed with high-pressure situations in mind including the signature red color and a breakaway lock system. Using a wide assortment of **Accessories**, these carts have the right product for any emergency situations.

Our QT Cart Emergency Package Model QT-EMERG-PKG includes the following accessories:

- VIV – IV Pole
- VDS – Defibrillator Shelf
- VCPR – Cardiac Board and Hooks
- VO2 – Oxygen Tank Holder
- VTR – Three-Sided Rail
- QT-AD3 – (2) 3” Drawer Divider Set
- AT-AD6 – (1) 6” Drawer Divider Set
- TOS-15 – 15’ Cord, 6 Outlets



MODEL # QT41B27K33669B  
Shown with Accessories



MODEL # QT41R27B33669R  
Shown with Accessories

## ANESTHESIA CARTS

Anesthesia Carts store all necessary tools in a single, compact portable location with a locking system designed to keep unauthorized personnel out. The blue color makes cart identification easy.

Our QT Cart Anesthesia Package Model QT-ANES-PKG includes the following accessories:

- VIV – IV Pole
- P-AB – Back Bridge
- VTB-1 – Two (2) Tilt Bins, 5-Bin System
- VGD – Glove Dispenser
- VWB – Waste Basket
- QT-AD3 – (2) 3” Drawer Divider Set
- AT-AD6 – (1) 6” Drawer Divider Set
- VTR – Three-Sided Rail
- TD-1 – Tape Dispenser



# SPECIALTY QT® CARTS & ACCESSORY PACKAGES



## ISOLATION CARTS

Isolation Carts are a convenient way to store PPE in an easy-access location. The bright yellow color is a constant reminder to use proper PPE and precautions for infection control when entering a patient room.

Our QT Isolation Package Model QT-ISO-PKG includes the following accessories:

- VWB – Waste Basket
- VGD – Glove Dispenser
- VTR – Three-Sided Rail
- VUH – (2) Utility Hooks with Accessory Rail
- VCR – Coat Rack

MODEL # QT44Y30K3999Y  
Shown with Accessories

## PEDIATRIC EMERGENCY CARTS

Pediatric Emergency Carts come with nine (9) drawers that are color code compatible with pediatric color-coding systems. A **Key Lock** and **Sliding Table Top** extension (max load of 17 pounds) are standard.

Our QT Pediatric Model QT44W30BPEDS2 includes the following:

- Lightweight Aluminum Design
- Aluminum Rails for Simple Accessory
- Overall Dimension – 20” D x 27” W x 43.5” H
- 30” Vertical Drawer Space
- Drawer Interior Dimensions – 15.25” D x 21.75” W
- Each Drawer Has Separate Breakaway Locks
- 5” Quiet Medical Casters – (1) Swivel, (2) Brake, (1) Swivel/Lock



MODEL #  
QT44W30BPEDS2

# SECURITY OPTIONS

## SECURITY OPTIONS

To help maintain patient safety and secure supplies, **QT Carts** come standard with **Key Locks**, however there are several lock options available. From **Breakaway Locks**, which offer quick access in critical situations, to **Proximity Keyless Entry** for the most secure applications, LogiQuip has the right lock system for every situation.

## AVAILABLE LOCK OPTIONS

- **Central Key Lock** - Comes standard on ALL carts.
- **Breakaway Lock** - Vertical bar for quick access during emergency situations.
- **Keyless Electronic Lock** - Touchpad for easy access.
- **HID Proximity or iClass Lock** - Touchless badge entry records access for increased security.




Central Lock  
Secures All Drawers

## ELOCK SPECS


- Whisper Quiet Operation
- 5,000 User Capacity
- Automatic Relock
- Key Override
- Dual Credential Capability
- Available WiFi Allows IT To:
  - Remotely Add/Delete Users
  - Automatically Download Audit Trails
  - Monitor Battery Level
- Requires LockView (VLS) For Remote Connection



eLock System With Proximity Switch  
(Optional)



MODEL # VLS  
Lock Management Software



HID (MODEL # VPC) and  
iClass Proximity Card







MODEL # VPT  
HID Proximity Tag

# CUSTOM CONFIGURATIONS & COLOR OPTIONS

## CUSTOM CART/CABINET CONFIGURATIONS

LogiQuip offers full customization of QT Carts from **Cabinet Size**, **Drawer Size**, and **Color**.


- Pick the **Cart Body** – Standard (27" wide) or Narrow (20" wide)
- Pick the **Body Color** –  White,  Red,  Yellow,  Blue
- Pick your **Drawer Size(s)** – 3", 6", 9"
- Pick your **Drawer Color** (see **Drawer Front Color Options** below)
- Add your **Accessories** (see pages 13-15)



Narrow QT Cart



Light Blue Drawer Color Option  
MODEL# QT38W24E3669LB



Standard QT Cart

## DRAWER FRONT COLOR OPTIONS

LogiQuip QT Cart **Drawer Fronts** are available in twelve (12) colors. **Body Colors** in White and limited quantities of Red, Blue and Yellow are available through the **Rapid Ship Program**. All **Drawer Colors** and **Sizes** are in stock and available through **Rapid Ship**.

(Actual color may vary, depending on your monitor setting or press run.)



Blue (B)



Cream (C)



Green (G)



Lavender (L)



Light Blue (LB)



Orange (O)



Pink (P)



Red (R)



Silver (S)



Taupe (T)



Yellow (Y)




White (W)


# LOGIQUIP® QT® CARTS

## QUIET STORAGE FOR PATIENT CARE

LogiQuip Quiet Treatment (QT) Carts comes standard with Key Lock for secure storage. With soft-closing drawers and hospital grade casters, our QT Carts help provide a quieter environment for patient care. Choose from our selection of Stock Accessories (pages 13-15) to customize your QT Cart.



Scan QR Code to see product details.



# STOCK ACCESSORIES (RAPID SHIP PROGRAM)

## STOCK ACCESSORIES

All Stock Accessories are available for immediate shipment through our Rapid Ship Program.



CPD



P-AB



QT-AD3, QTN-AD3



QT-AD6, QTN-AD6



TD-1



VAB

### RAPID SHIP


#### STOCK ACCESSORIES

MODEL	DESCRIPTION
CPD	Computer Drawer
P-AB	Back Bridge
QT-AD3	3" Standard Cart – 3 short, 6 adjustable, 2-3" strips, 6 T-clips
QT-AD6	6" Standard Cart Set – 3 short, 6 adjustable, 2-6" strips, 4 T-clips
QTN-AD3	3" Narrow Cart Set – 3 short, 6 adjustable, 2-3" strips, 6 T-clips
QTN-AD6	6" Narrow Cart Set – 3 short, 6 adjustable, 2-6" strips, 4 T-clips
TD-1	Tape Dispenser
TOS-15	15' Cord, 6 Outlets
VAB	Accessory Bar
VAD-1	Locking Medicine Drawer
VBS-1†	Small Fixed Bins and Accessory Bar Each Bin 7" x 4" x 3"

†Requires Back Bridge (P-AB)



VAD-1



VBS-1†



# STOCK ACCESSORIES (RAPID SHIP PROGRAM)



VCPR



VCR



VDS



VFH



VGB



VGD



VHL



VGD-3




VIV


## RAPID SHIP

### STOCK ACCESSORIES

MODEL	DESCRIPTION
VCPR	Cardiac Board and (2) Mounting Hooks
VCR	Coat Rack
VDS	Defibrillator Shelf
VFH	File Holder and Accessory Bar 6" x 15" x 22"
VGB	Garbage Can With Lid, Basket and Accessory Bar
VGD	Gloves Box Holder and Accessory Bar
VGD-3	Triple Glove Box Holder, Accessory Bar (2)
VHL	Push Handle with Accessory Bar
VIV	IV Pole
VMS	Multi-Storage Basket and Accessory Bar 7" x 13" x 5"
VMT-1	Medicine 15" x 22" x 2" with 9 compartments



VMS



VMT-1

# STOCK ACCESSORIES (RAPID SHIP PROGRAM)



VMT-2




VO2



VSC



VSH



VTR



VTB-1†

VTB-2†



VTS†



VUH



VWB



VWC

**RAPID SHIP**

STOCK ACCESSORIES	
MODEL	DESCRIPTION
VMT-2	Medicine 15" x 22" x 4.25" with 9 compartments
VO2	O2 Holder and Accessory Bar
VSC	Sharp Container Set, Basket and Accessory Bar
VSH	Suction Machine and Accessory Bar
VTB-1†	Large Tilt Bins (5)
VTB-2†	Medium Tilt Bins (12)
VTR	Three Sided Rail
VTS†	Top Shelf and Accessory Bar for QT Carts
VUH	Utility Hook Rack (2 Hooks) and Accessory Bar
VWB	Waste Basket and Accessory Bar
VWC	Garbage Can With Rotating Lid, Basket and Accessory Bar

†Requires Back Bridge (P-AB)

