

Intelli-Flip™

INVENTORY MANAGEMENT



INTELLI-FLIP™ SYSTEM

“ ...the LogiQuip (Intelli-Flip) system is extremely simple and we feel will be a big benefit to our team as well as Nursing... feedback and benefits of this system were very positive from Nursing and Supply Chain staff... my director is looking to push the system out to all sites in the near future. ”

Ben McGlashen, Regional Manager Supply Chain Site Operation



INTELLI-FLIP™ SYSTEM

Every day, Healthcare Professionals are asked to **Do More With Less**. With **Intelli-Flip**, you can do just that. LogiQuip's patented **Intelli-Flip System** is a simple, intuitive tool that automates inventory management, **saving you up to 50 minutes per day per PAR location in supply management tasks**, helping to ensure you have the supplies you need when you need them. **No scanning. No weight calibration. No wires to work around.**

SIMPLE & RESPONSIVE

An **Intelli-Flip Monitor** is attached to each **Kanban Flipper**. When the caregiver has used all of the supplies in the first bin, they simply lower the **Kanban Flipper** to access supplies in the second bin, and a signal is sent to the **Portal** that inventory needs to be replenished. **No scanning required.**

INTELLI-FLIP PORTAL

Track inventory, velocity trends, days on hand, stock levels, and system status updates on Monitor function through a secure **Portal** on Intelli-Flip.net. **No need to go to the PAR location.**

ERP INTEGRATION

Easily create orders with **ERP Integration**. **Intelli-Flip** can be integrated into your **ERP System** to automate ordering, reducing stock-outs and expired products.

ACCESS CONTROL

Users will have access to system data and adjustments based on assigned level of responsibility. Training by user level will be done upon install.

INTELLI-FLIP CONNECT

Intelli-Flip Connect is an exclusive feature that allows the system to move to a cellular connection if the Wi-Fi signal is lost. This helps to ensure you never lose information.

IMPLEMENTATION

LogiQuip has teams of experienced professionals that will install and set-up the **Intelli-Flip System**, allowing you to experience the benefits of the system quickly!

Intelli-Flip™ is Simple & Responsive

— REDUCES SUPPLY CHAIN COSTS AND STOCKOUTS —

NO

SCANNING
WEIGHT CALIBRATION
WIRES

SAVE 50 MINUTES PER PAR LOCATION PER DAY* ON SUPPLY MANAGEMENT

ERP INTEGRATION

NO NEED TO GO TO THE
PAR LOCATION

*Assumes 272 SKUs per PAR location with an average time savings of 11 seconds per SKU based on a case study done with Intelli-Flip.

INTELLI-FLIP™ MANAGES SUPPLY INVENTORY

Managing supply inventory has never been easier! **Intelli-Flip** replaces the traditional two-bin system and **automates supply tracking**. A **Monitor**, attached to the **Kanban Flipper**, placed in the center of the bin signals that PAR levels are low when it is pulled down to access supplies behind the **Kanban Flipper**. Inventory can be easily viewed and managed through a secure Portal on Intelli-Flip.net.

FULLY STOCKED SHELVES

When shelves are fully stocked, all of the **Kanban Flippers** are in the **UP** position.



GET SUPPLIES

Easily obtain the supplies from the first bin when you need them.



PULL DOWN THE FLIPPER WITH INTELLI-FLIP

Once the supplies in the first bin are gone, simply pull down the **Kanban Flipper** and a signal will be sent to the **Portal** to indicate that more stock is needed.



CONTINUE TO GET SUPPLIES

Continue to use supplies from the second bin, behind the **Kanban Flipper**, until the bin is restocked.



INTELLI-FLIP™ PORTAL

Using the secure **Intelli-Flip Portal** on Intelli-Flip.net allows you to manage and requisition supplies without scanning or going to PAR Locations. **Intelli-Flip** gives you access that is integrated into your **ERP System** for easy ordering.

MANAGE INVENTORY

Fully automate supply tracking with **Intelli-Flip**. The **Intelli-Flip Portal** tells you which supplies need to be replenished at each PAR Location without scanning or visiting the PAR locations.

ERP INTEGRATION


Easily create requisition files to upload into your **ERP System** for order generation.

TRACK PRODUCT ACTIVITY TRENDS

Track **Product Activity Trends** and order accordingly based on the data from your **Portal**.

MONITOR SUPPLIES

Monitor which supplies are used most, how many days on hand of supplies, when stock levels dip, and when supplies exceed expected velocity so that you can react accordingly.



The screenshot shows the 'Details' page of the Intelli-Flip Portal, displaying a table of supply items. The table includes columns for Status, Name, Description, SKU, Timestamp, Mapping, Birt., Flipper, Qty, and UoM. The items listed include Odor Eliminator, Skin Macker, Facial Tissue, Irrigation Tray, Denture Cups, Enalapril Collection Set, Ice Bags, Thermometer, Mucus Trap, and Suction Collector.

Status	Name	Description	SKU	Timestamp	Mapping	Birt.	Flipper	Qty	UoM
Good	Odor Eliminator	DEODORANT ODT 10Z UNSCENT PMP 87L SIMPLY FRESH CARRINGTON	135639	1/26/2024 10:02 AM	Aisle 1 Shelf 3	Rack A Bin 8	Good Down	12	Each
Good	Skin Macker	MARKER SUR RUL AND BLNK LBL 6IN SKIN REG TIP W/ L RSK RESVR	187508	1/26/2024 9:48 AM	Aisle 1 Shelf 3	Rack A Bin 9	Good Up	6	Each
Good	Facial Tissue	TISSUE FACE WS.76L7IN STD 2 PLY RECYCL PAPER	175329	1/26/2024 9:57 AM	Aisle 1 Shelf 4	Rack A Bin 1	Good Up	27	Box
Good	Irrigation Tray	TRAY IRIG 60ML CATH TIP PIST SVR THMB RING ELONG ALC PD SQ	104133	1/26/2024 9:53 AM	Aisle 1 Shelf 4	Rack A Bin 10	Good Up	6	Each
Good	Denture Cups	CLIP ASP MUSCDD01 DENTURE 8 OZ NEON GREEN	233627	1/26/2024 9:59 AM	Aisle 1 Shelf 4	Rack A Bin 2	Good Up	13	Each
Good	Enalapril Collection Set	KIT COLLECTION BD EDWAB REGULAR FLOOR 39948	108513	1/26/2024 10:00 AM	Aisle 1 Shelf 4	Rack A Bin 3	Good Up	12	Each
Good	Ice Bags	BAG THER W/50L CATH RECTANG STD ICE REFILLABLE W/ CLMP CLR ADJ E STRP MULT USE VELC STRP	187045	1/26/2024 10:01 AM	Aisle 1 Shelf 4	Rack A Bin 4	Good Up	3	Bag
Good	Thermometer	THERMOMETER DGT ORAL STD FC 30 2ND 1 PC SHMP CA 1447CA	139203	1/26/2024 10:01 AM	Aisle 1 Shelf 4	Rack A Bin 5	Good Up	5	Each
Good	Mucus Trap	TRAP MUCUS PLASTIC SCREW CAP TUBE ID LABEL STERILE LATEX FREE CLEAR	187261	1/26/2024 10:03 AM	Aisle 1 Shelf 4	Rack A Bin 6	Good Up	4	Each
Good	Suction Collector	SYSTEM SPUTUM COLLECTION FUNNEL HINGE UID CONICAL TUBE REMOVABLE BASE FALCON 50 ML 45868	105592	1/26/2024 9:49 AM	Aisle 1 Shelf 4	Rack A Bin 7	Good Up	3	Each

INTELLI-FLIP™ PRODUCT INTEGRATION

Intelli-Flip integrates with our existing patented **ParWall®**, **ParWire®**, **LQS Wire Shelving**, and **Kanban Flipper** products for easy installation. Once **Intelli-Flip** is installed, inventory levels for all supply locations using **Intelli-Flip** can be viewed and managed from a secure **Portal** at Intelli-Flip.net.

PARWALL

LogiQuip's **ParWall Systems** are designed to maximize storage space, provide high product visibility, and increase usable space. **ParWall Systems** can be **Freestanding (Stationary, Mobile)** or fixed to a wall (**Lean Wall®** or **Louvered Panel**). These **Frames** can accommodate eight (8) **Baskets** vertically and accept **Slanted Baskets** for easy access to supplies.



PARWIRE

The LogiQuip **ParWire System** with our **Fine Mesh Basket** design helps to increase storage capacity by up to 40%. These strong and durable shelves come in a variety of sizes and are compatible with LogiQuip's 1-inch round **Posts** for easy customization. This design is available in **Standard, Slanted, or Glide** options. Our **ParWire Fine Mesh Glide Basket** extends out 50% for full visibility and easy access to stored items.



LQS WIRE SHELVING

In addition to saving space, our **LQS Wire Shelving** and **Posts** begin with a **Carbon Steel Core** that is **Double Chrome Plated** for strength and durability that lasts. With a wide variety of configurations available to meet any of your organizational needs, our **LQS Wire Shelving Line** is flexible enough to be used virtually anywhere within your facility for easy and efficient storage.



KANBAN FLIPPERS

The **Intelli-Flip Monitor** easily attaches to our **Kanban Flipper**, and is triggered by the natural act of accessing the stock from the second bin.



INTELLI-FLIP™ SYSTEM INSTALLATION

Let us do the heavy lifting! From your first consult and throughout your day-to-day experience with **Intelli-Flip**, LogiQuip's **Services Team** is there. Our experienced staff will work with you to review your existing layout and PAR levels, implement the final solution, and train your staff, making this a turnkey solution that will help ensure you have the supplies you need when you need them.

REVIEW EXISTING LAYOUT AND PAR LEVELS

The LogiQuip **Services Team** is ready to work with you to review your existing space, supply usage, and your future goals to determine the best design for your unique situation.

- Obtain PAR Lists
- Review Floor Plans
- Cart Layouts and Photos
- Work with Staff to Determine Issues to Address

IMPLEMENT NEW DESIGN

Our **Services Team** is dedicated to providing on-time delivery and efficient installation of your **Shelving units** and **Intelli-Flip System** to get you up and running as quickly as possible.


- Build-Out Shelving Units
- Install **Intelli-Flip Monitors** and Assign to Products
- Install Data Consolidator and Connect to [Intelli-Flip.net](https://intelli-flip.net)

TRAIN STAFF

Start realizing the benefits of the **Intelli-Flip System** from day one with full system training with a member of our **Services Team**.

- Train System Administrator
- Train Department Staff
- Train Nursing Staff





www.logiquip.com

(800) 665-3760 | customerexperience@logiquip.net

5.2024