



Case

CARTS



WHY LOGIQUIP®?



Healthcare exclusive since 1992, LogiQuip's **Inventory Management Solutions** are designed for Healthcare Professionals looking for a better way to manage supply chain tasks. With **Expert Advice**, **Innovative Solutions**, and **360° Service** that doesn't stop after the sale, LogiQuip's sole purpose is to ensure Healthcare Professionals have the supplies they need when they need them. Partnering with you to design and implement the right solutions, we give you time back to focus on what's important, Patient Care.

EXPERT ADVICE

With an in-depth understanding of the industry and a focus on **Infection Control** and **Lean Processes**, we partner with you to ensure we design the right solutions for you.

INNOVATIVE SOLUTIONS

Designed with an eye on the horizon, LogiQuip has **Innovative Solutions** that evolve with your facility. Our flexible and efficient options help to ensure you have the supplies you need when you need them.

360° SERVICE

LogiQuip's experienced **Sales Representatives** and **Award-Winning Customer Support Team** provide unparalleled service before, during, and after the sale.



LogiQuip is committed to fulfilling customer orders as quickly as possible which is why we have hundreds of products in stock, ready to ship every day. Please contact LogiQuip's Customer Experience Team or your Sales Representative to learn more.

WITH LABOR SHORTAGES AND THE EVER-INCREASING DEMAND TO DO MORE WITH LESS, the **LogiQuip Services Team** can help to decrease that burden with turnkey solutions to help ensure you have the supplies you need when you need them.

PROJECT MANAGEMENT

From our initial meeting with you through implementation, a LogiQuip **Services Team Member** will be assigned as your **Project Manager**, working with you every step of the way. Your **Project Manager** will serve as your single point of contact during the project life cycle. Together we will develop and implement project milestones until the work is complete. Contact your **Sales Representative** or call **(800) 665-3760** for more information.

ON-SITE ANALYSIS

Working with your LogiQuip **Representative**, an on-site analysis will be done to ensure that the **Project Manager** has a complete understanding of your current environment. During our on-site analysis we will review your current state, identify gaps, and align on the goals of the project.

PLANNING

Using a mix of innovative products, Lean principles, tried and true methods, and out of the box thinking, LogiQuip will provide a holistic plan with turnkey solutions to minimize downtime while decreasing your storage footprint and increasing PAR levels.

INSTALLATION & CONVERSION

LogiQuip is dedicated to providing on-time delivery, efficient installation, and supply conversion to minimize down time and time on site to get you up and running as quickly as possible.

ONGOING SUPPORT

With the project complete, we want to help you sustain the changes that have been made. With **Annual Inventory Reviews**, **Preventative Maintenance Plans**, and **Field Service** for in warranty and out of warranty repairs, we have you covered to keep things running smoothly.

LOGIQUIP® CASE CARTS

Designed specifically to protect sterile supplies as they move from central processing to surgical suites, LogiQuip **Case Carts** provide a safe, durable, and efficient means of transportation. All **Stainless Steel Bolted** and **Over-the-Road Case Carts** have a full wrap-around **Bumper**, **6" Casters** for easy maneuverability, and accept height adjustable **Shelving** (**Shelves** not included).

STAINLESS STEEL CONSTRUCTION

An easy-to-clean, corrosion-resistant **Stainless Steel** body on all **Closed** and **Over-the-Road Case Carts** provides strength and durability in the harshest environment.

6-INCH CASTERS

These durable **Casters** can withstand high temperatures and high pressure when going through an automated washer. Sealed ball bearings allow these **Casters** to quietly move through your facility.

ACCESSORIES & COMPONENTS

With a variety of bolt-on accessories and components including **Card Holders**, **Clean/Dirty Indicator**, **Handles**, and **Bumpers**, you can customize these **Case Carts** to help streamline processes.

FULLY SERVICEABLE CONSTRUCTION

Case Carts are abused daily, which can lead to damage. With **Bolted** (vs. welded) **Construction**, LogiQuip **Case Carts** are easily serviced to reduce downtime and extend the life of your **Case Carts**.



BOLTED CASE CARTS



OVER-THE-ROAD CASE CART

LIMITED 5-YEAR WARRANTY

LogiQuip products are covered by a **Limited 5-Year Warranty** from date of receipt against defect in materials or workmanship under normal use and service, excluding **Electronic Components**, **Cart Covers**, **Custom Products**, **LogiCell Tambour Doors**, **Label Holders**, and **Casters**. During this period, LogiQuip will, in its sole discretion, repair defects in materials or workmanship that existed when the product was received or replace the product. Visit [LogiQuip.com/Warranty](https://www.logiquip.com/Warranty) for more details.

LOGIQUIP BOLTED STAINLESS STEEL CASE CARTS

BOLTED STAINLESS STEEL CASE CARTS

LogiQuip Bolted Case Carts are easily serviceable, reducing downtime and extending the life of your Case Carts. Large enough to accommodate oversized instrument sets, these Case Carts come standard with Ergonomic Handles, 270° Door Swings, 6" Casters, and Wrap Around Bumper. The Stainless Steel Cart Body is easy-to-clean, has no water retention points, and comes with a no drip "tilted" pan base. Shelving options are sold separately.

BOLTED STAINLESS STEEL CASE CARTS		
MODEL	EXTERIOR D x W x H - IN.	INTERIOR D x W x H - IN.
LQCC3630	27 x 36 x 30	25 x 30.5 x 22
LQCC3640	27 x 36 x 40	25 x 30.5 x 32
LQCC3655	27 x 36 x 55	25 x 30.5 x 47
LQCC4640	27 x 46 x 40	25 x 40.25 x 32
LQCC4655	27 x 46 x 55	25 x 40.25 x 47
LQCC5255	28 x 52 x 55	25 x 46.5 x 47




MODEL # LQCC3640

SHELF OPTIONS FOR BOLTED CASE CARTS

There are three (3) Roll-Out Shelving options available for Bolted Stainless Steel Case Carts; Perforated Stainless Steel Shelves, Solid Stainless Steel Shelves, and Wire Case Shelves. All Shelves are sold separately and come with Adjustable Rails. Maximum weight capacity is 200 pounds for Stainless Steel Shelves and 150 pounds for the Wire Shelves.

BOLTED CASE CART SHELF OPTIONS		
MODEL	DESCRIPTION	D x W - IN.
LQCCS36PR	Roll-Out Perforated Stainless Steel Case Cart Shelf for LQCC3630, LQCC3640, and LQCC3655	23 x 26
LQCCS36SR	Roll-Out Solid Stainless Steel Case Cart Shelf for LQCC3630, LQCC3640, and LQCC3655	23 x 26
LQCCS36WR	Roll-Out Wire Case Cart Shelf for LQCC3630, LQCC3640, and LQCC3655	23 x 26
LQCCS46PR	Roll-Out Perforated Stainless Steel Case Cart Shelf for LQCC4640 and LQCC4655	23 x 36
LQCCS46SR	Roll-Out Solid Stainless Steel Case Cart Shelf for LQCC4640 and LQCC465	23 x 36
LQCCS46WR	Roll-Out Wire Case Cart Shelf for LQCC4640 and LQCC465	23 x 36
LQCCS52PR	Roll-Out Perforated Stainless Steel Case Cart Shelf for LQCC5240 and LQCC5255	23 x 42
LQCCS52WR	Roll-Out Wire Case Cart Shelf for LQCC5240 and LQCC5255	23 x 42



Roll-Out Perforated Stainless Steel Shelf



Roll-Out Solid Stainless Steel Shelf

Roll-Out Wire Shelf (not pictured)

OVER-THE-ROAD CASE CARTS

OVER-THE-ROAD CASE CARTS

Building on our **Standard Case Carts**, LogiQuip's **Over-The-Road Case Carts** provide a rugged, durable and safe way to transport sterile supplies when Central Processing is off-site. All welded, stainless construction, heavy-duty 6" x 2" washable casters and optional accessories, help to create a clean and safe environment.

All Over-The-Road Case Carts are available as a Custom Build.

- Exterior Dimension: 30" D x 43.75" W x 57" H (8.28 sq. ft.)
- Shelf Size: 26" D x 39" W
- Capacity: Six (6) full-size containers
- Container Capacity: Available for 12, 18, and 24 containers
- Unique container retention system
- Storage Space: Top 5"H; Bottom 14.5" H
- Casters: Two (2) directional locking; Two (2) swivel with brake
- Bumper System: Full wrap around on top and bottom of cart
- Door System: Overlapping, self-sealing with 95° cart seal and double locking key
- Security Tab: Breakaway



MODEL # OTR0100-01-KT Easy slide CLEAN/DIRTY Indicator eliminates guesswork and errors



MODEL# CC-FLAG3 Flag System

OVER-THE-ROAD CASE CART CASTERS/ACCESSORIES	
MODEL	DESCRIPTION
OTR0100-01-KT	Clean/Dirty Indicator
CC-FLAG3	Flag System

CASE CARTS ACCESSORIES & CASTERS

CLOSED CASE CART ACCESSORIES

Accessories for LogiQuip's Closed Case Carts are shown at right and on page 6.

CLOSED CASE CART ACCESSORIES	
MODEL	DESCRIPTION
BO-0204	Pick List Card Holder - 8.5" x 11" Stainless
CC-FLAG3	Flag System (shown on page 6)



MODEL# BO-0204
Card Holder

CLOSED CASE CART CASTERS

Our Standard Case Cart Casters are 6" x 1.5" stainless steel, all swivel (one is directional locking and two are braking), and have the following:

- Sealed bearings - upper and lower
- Grease fitting - bearings and axle
- High temp wheels - poly
- Cart washable
- Optional sizes available upon request



MODEL#
CCCAST6




MODEL#
CCCAST6B



MODEL#
CCCAST6DL

CLOSED CASE CART CASTERS	
MODEL	DESCRIPTION
CCCAST6	Cart Caster - Swivel
CCCAST6B	Cart Caster - Swivel with Brake
CCCAST6DL	Cart Caster - Directional Lock





www.logiquip.com

(800) 665-3760 | customerexperience@logiquip.net

LCCC-04-26